

Product Information

1.	Article number:	9051015		
2.	Product name:	Askina® Brauncel		
3.	Product description:	Cellulose swabs of high-bleached surgical cellulose		
4.	List of materials:	see attachments		
5.	Form of delivery:	☐ sterile	☒ non-sterile	☐ non-sterile/bulk
6.	Information about sterilization:	see attachments		
7.	Latex (natural rubber):	☐ yes	☒ no	
8.	DEHP (according to REACH regulation)	☐ yes	☒ no	
9.	Shelf life:	☒ 5 years	☐ others	☐ not applicable
10.	Instructions for use:	☐ yes	☒ not applicable	
11.	Bioburden:	our non-sterile products do not exceed the Bioburden value of 200 CFU (colony forming unit)		
12.	Is the product CE marked?	☒ yes	☐ no	
13.	Identification number of the Notified Body	not applicable		
14.	Classification according Council Directive 93/42/EEC concerning Medical Devices Annex IX:			
		class: I	rule: 4	
15.	Attachments, if applicable:	ISO 13485 Certificate Declaration of Conformity (DoC) EC Certificate List of materials Confirmation Sterilization / Re-Sterilization Instructions for use (IFU)		

All data refer to typical single values and may be subject to alterations, valid without signature.

Edition: September 2018