Water bath protective media Aqua Stabil 6 x 100 ML



JULABO AQUA STABIL is a water disinfectant for use in open baths. This disinfectant prevents contamination of the bath fluid and ensures a longer service life when performing temperature control applications with water. The germicidal effect prevents the formation of algae, bacteria, and microorganisms.

Benefits:

- Excellent heat conductivity
- Minimum odor

Application: Add 2 ml of AQUA STABIL per liter of water. Replace water and AQUA STABIL if the blue color fades, or after three months at the latest. For the expiration date, see label on reverse side.

Order number	**8940006**