Allgemeine Produktdaten

ratiolab® Kryo-Boxen mit Rastereinsatz
Autoklavierbar, temperaturbeständig von -90°C bis +121°C. 81 Stellplätze im 9 x 9 Raster nummerisch codiert. Stabiler Scharnierdeckel mit sicherem Schnappverschluss. Stapelecken für den sicheren Transport. Ablauföffnungen für Tauflüssigkeit am Boden. Herausnehmbare Rastereinsätze. Die Kryo-Boxen sind ohne Rastereinsatz auch für größere Gefäße verwendbar.

General Product Information

ratiolab® Cryo boxes with Grids
Autoclavable, temperature-resistant from -90°C to +121°C. 81 places, 9 x 9 grid, numerically coded. Robust hinge with safe snap-on lid. Stackable for safe transport. Drain bores at the bottom for dew liquids. Removable grid inserts. Cryo Boxes can also be used for larger vessels by removing the grid inserts.

Material / *Material:* Polypropylen / *Polypropylene*

Warentarif-Nummer / *Intrastat Number- HS Code:* 39269097

CE-Kennzeichnung / *CE-marking:* Ja/*Yes* ☐ Nein/*No* ☒

Sterilisiert / *pre-sterilized*: Ja/*Yes* ☐ Nein/*No* ☒

Farbe *Color*	B x T x H mm *W x D x H mm*	für Röhrchen Ø mm For Tube Ø mm	Verpackungseinheit *Quantity per Pack*	Raster *Grid*	Bestell-Nr. *Order No*
für Kryo - Röhrchen 1,2 - 2 ml / *for Cryo Tubes 1.2 - 2 ml*					
natur / *natural*	133 x 133 x 52	12,8	5	9 x 9	51 20 010
gelb / *yellow*	133 x 133 x 52	12,8	5	9 x 9	51 20 011
rot / *red*	133 x 133 x 52	12,8	5	9 x 9	51 20 012
grün / *green*	133 x 133 x 52	12,8	5	9 x 9	51 20 013
blau / *blue*	133 x 133 x 52	12,8	5	9 x 9	51 20 014
schwarz / *black*	133 x 133 x 52	12,8	5	9 x 9	51 20 015

Farbe *Color*	B x T x H mm *W x D x H mm*	für Röhrchen Ø mm *For Tube Ø mm*	Verpackungseinheit *Quantity per Pack*	Raster *Grid*	Bestell-Nr. *Order No*
für Kryo - Röhrchen 3 - 4 ml / *for Cryo Tubes 3 - 4 ml*					
natur / *natural*	133 x 133 x 75	12,8	5	9 x 9	51 20 020
gelb / *yellow*	133 x 133 x 75	12,8	5	9 x 9	51 20 021
rot / *red*	133 x 133 x 75	12,8	5	9 x 9	51 20 022
grün / *green*	133 x 133 x 75	12,8	5	9 x 9	51 20 023
blau / *blue*	133 x 133 x 75	12,8	5	9 x 9	51 20 024
schwarz / *black*	133 x 133 x 75	12,8	5	9 x 9	51 20 025

erstellt: Anna Chrupala
06.08.2019 A. Chrupala Praktikantin
Datum, Unterschrift / *Date, Signature*

geprüft: E. Tagik
06.08.2019 E. Tagik QS/ Technik
Datum, Unterschrift / *Date, Signature*