Technical sheet

May-Grünwald Solution

Code product: **320070**

Synonym:

Methylene Blue-Eosine Solution

Specifications:

- **Appearance:** Clear dark blue Solution
- **Absorption Maximum in Methanol:** 522 nm ≤ λ_1 ≤ 526 nm
  651 nm ≤ λ_2 ≤ 653 nm
- **pH (50/50 in Demineralised Water):** 6.5 -7.5

Applications:

Haematology: used in May-Grünwald Giemsa technique for differential staining of medullary and blood smears.

Warning:

This reagent must be transported, stored and used at a temperature between 15 and 25 °C

Updated: 25/04/2013

RAL Diagnostics - Site Montesquieu - 33651 MARTILLAC Cedex – France - Tel : 33+557 960 404 - Fax : 33+557 960 405 - www.ral-diagnostics.fr