
:\DWGS\CSD-W216909


THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF DWK LIFE SCIENCES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES INC. ALL RIGHTS RESERVED.©
ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

THIRD ANGLE PROJECTION

PRODUCT COMPONENTS

ITEM	QTY	DESCRIPTION





FINISH DETAIL
SCALE 2 X 1

JAR MATERIAL: TYPE III AMBER SODA–LIME GLASS
CAP MATERIAL: WHITE POLYPROPYLENE
C LINER MATERIAL: 1.02 THICK WHITE POLYETHYLENE FOAM WITH 0.13 NATURAL PTFE–FACED LAMINATED TOGETHER

3. WEIGHT: BOTTLE: 113.4g (APPROX.) | CAP: 5.5 ±1.5g. (EXCLUDES LINER)
2. DEPRESSED THREADS AT SEAMS.
1. DIMENSIONS ARE IN MILLIMETERS.
----NOTES----

TITLE: 4oz JAR, STR, TYPE III CLR, PP/PTFE LNR B
CAPACITY: 4.62 ±.078oz (136.64 ±2.3cc) OVERFLOW
FINISH: G.P.I. 58–400–405 U'CUT B
SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.

DATE DRAWN 02/12/07 | PROJ # | – | REVISED TO DWK STDS | 06/15/19 EL
DRAWN BY B. LLOYD | APPROX. WGT. SEE NOTE 3 | C | 2014051410392381 | 06/05/14 TP
CHECKED BY TRJ 2/13/07 | AUTHORIZED CJM 2/15/07 | REV | E.C.N. # | DATE

DURAN WHEATON KIMBLE
1501 North 10th Street | Millville, NJ 08332–2038 | USA

DWG. NO: B– CSD–W216909