


29.11.2022

PRODUCT SPECIFICATIONS SHEET

Mortar with lip unglazed

Jizerská porcelánka s.r.o. certifies that the specified products meet all specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

Sign	JIPO Code	D (mm)	H (mm)	V (ml)	
211a/0	641 331 211 000	51	32	25	
211a/1	641 331 211 100	56	36	30	
211a/1/0	641 331 211 110	63	41	50	
211a/1a	641 331 211 120	76	45	75	
221a/2	641 331 211 200	83	50	110	
211a/2a	641 331 211 220	90	56	160	DIN 12 906
211a/3	641 331 211 300	105	64	220	
211a/3a	641 331 211 320	125	71	400	DIN 12 906
211a/4	641 331 211 400	139	82	500	
211a/5	641 331 211 500	150	90	650	DIN 12 906
211a/6	641 331 211 600	180	104	1000	DIN 12 906
211a/7	641 331 211 700	224	127	2250	
211a/8	641 331 211 800	268	160	4300	



Jizerská porcelánka s.r.o
Údolní 138
Desná v Jizerských horách

VAT : CZ 46710663
e-mail : info@jipo-desna.cz
http://www.jipo-desna.cz

tel : 00420 483 369 710
fax : 00420 483 383 327



29.11.2022 **PRODUCT SPECIFICATIONS SHEET**

Mortar with lip glazed

Intended uses: Laboratory applications

Material: Porcelain C110

Shape: Deep form with spout and without graduation

Compliant to: ISO material for the manufacture lies within the group C110 according to ČSN EN 60672 –3
Product dimension tolerances are controlled by DIN 40 680
Acid resistence acc. to ISO 1775

Properties: Unglazed products can be used at temperatures of up to 1350°C / 2462°F and glazed ones at temperatures of up to 1000°C / 1,832°F (Δ t 230°C)
Products have high resistance to acids and bases
Density 2.4 g/cm^3
Dending strength of unglazed products (70MPa)
Thermal conductivity (W.m^{-1}.K^{-1}) = 1.3

The compliance of demands is proved by regular tests and followed noticed tests from an independent laboratory:

Test according ČSN EN 60672-2
Test resistance to acids and bases
Complies with Council Directive 84/500 EEC
Complies with Regulation No. 1935/201 of the European Parliament and of the Council

Disclaimer: *Jizerská porcelánka s.r.o. wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use. If there's any doubt ask us.*

Jizerská porcelánka s.r.o
Údolní 138
Desná v Jizerských horách

VAT : CZ 46710663
e-mail : info@jipo-desna.cz
http://www.jipo-desna.cz

tel : 00420 483 369 710
fax : 00420 483 383 327