
CAT. NO.....:	155-0020/-0045		157-0020/-0045		156-4020/-4045		158-0020/-0045	
PORE SIZE (um).....:	0.2/0.45		0.2/0.45		0.2/0.45		0.2/0.45	
UPPER CAP. (mL).....:	150		250		500		500	
RECEIVER CAP. (mL).:	150		250		500		1000	
MEMBRANE COLOR.....:	WHITE		WHITE		WHITE		WHITE	
MEMBRANE DIA (mm)..:	50		50		75		75	
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A (UPPER I.D.).....:	3.15	80.0	3.15	80.0	3.80	96.5	3.80	96.5
B HEIGHT.....:	6.15	156.2	7.9	200	9.8	249	11.4	289
C DIA (BASE).....:	3.50	88.9	3.50	88.9	3.86	98.0	4.48	113.8
D OVERALL.....:	3.7	94	3.7	94	4.0	102	4.2	106
NOM. WT. (gms).....:	78		95		155		205	

MATERIALS	
UPPER CONTAINER	PS
COVER	PS
COLLAR	YELLOW PS
RECEIVER	PS
STORAGE CLOSURE	WHITE HDPE
ADAPTOR	WHITE LDPE
GASKET	FOAMED LDPE
MEMBRANE	SFCA
VENT PLUG	CA
PRE-FILTER	GLASS FIBER (156-XXXX & 158-XXXX ONLY)



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JGJ/MG	APP'D	PVC 6/25/99	PSS 00-15X-XXXX-1	REV. 0699
CK'D	PVC				

DESCRIPTION:

TISSUE CULTURE FILTER UNITS