PRODUCT DATA SHEET

STERICAN 0,80X80 G21X3 1/8 INTRAMUSCULAR (4665465)

Single-use hypodermic needles



Attribute Name	Attribute Value
Brand name	Sterican®
Product name (HC)	Sterican®
Product type	Needle
Application	Injection, Blood sampling
Connection	Luer-Lock/female
Type of bevel	long
Outer-ø cannula	0.80 mm
Gauge (cannula)	21
Color	green
Length of cannula (metric)	80.00 mm

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com

Attribute Name	Attribute Value
Length of cannula (inch)	3 1/8"
Residual volume	0.074 ml
Safety clip	No
Type of filter	No filter
Type of packaging	Blister
Product (unpack.) contains PVC	No
Product (unpack.) cont. DEHP	No
Product (unpack.) cont. Latex	No
Packaging contains PVC	No
Packaging contains DEHP	No
Packaging contains Latex	No
Product contains PVC	No
Product contains DEHP	No
Product contains Latex	No
Product contains Nickel	Yes
Sterility	sterile
Type of Sterilization	EO
Disposable / Reusable	disposable