Diethyl Ether, Anhydrous

Baker

Stabilized with BHT

Material No: 8106.9025
Effective from: 2013/07/22

Product Specification

Meets Reagent Specifications for testing Ph. Eur. monographs

Test	Specification
Acidity	passes test
Aldehydes	passes test
Distillation Range	34.0-35.0°C
Non-volatile Matter	<= 20 mg/l
Peroxide	passes test
Relative Density	0.714-0.716
Substances with a Foreign Odour	passes test
Water (H2O)	<= 2 g/l

For Laboratory, Research or Manufacturing use.

Storage conditions:

Store in an airtight container.
Store protected from light.
Store at a temperature of 8°C to 15°C.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000