Code: 44330000, 44300000 y 44310000
Wash bottles

Catalogue page 2014: 249

Code	Description	Case Qty	Case Weight. (kg)	Case Vol. (m^3)
44330000	Wash bottle 250ml	90	5.30	0.078
44300000	Wash bottle 500ml	60	4.85	0.080
44310000	Wash bottle 1.000ml	45	5.12	0.106

DESCRIPCIÓN

Made of a combination of high and low density polyethylene (HDPE and LDPE).
Transluscent bottle with HDPE cap (blue), delivery tube and seal cap (natural).
Bottle caps are ribbed for a better handling.
High resistance to all liquids, acids and bases.
Include a large neck for safe and easy filling.



MEASURES

Code	Weight g	Body Ø mm	Body Height mm	Neck Height mm	Thread Ø	Interior Ø
44330000	32 (±1.60)	66.80 (±0.80)	157.00 (±1.50)	30.40 (±0.40)	36.50 (±0.30)	19.80 (±0.20)
44300000	46 (±2.30)	78.50 (±0.80)	176.40 (±1.50)	30.40 (±0.40)	36.50 (±0.30)	19.80 (±0.20)
44310000	72 (±3.60)	95.50 (±0.80)	211.00(±1.50)	30.40 (±0.40)	36.50 (±0.30)	19.80 (±0.20)

Code	Total volume ml
44330000	410 (±10)
44300000	590 (±10)
44310000	1.070 (±10)

Elaborated	Verified	Authorised
Gemma Riba Marketing Department	Anna Mir Technical, Development and Quality Director	Monica Torras Sales and Marketing Director