
CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
002015	1	2	4	8	16	32	64	PC	PP*

* 53B CLOSURE - 1000 SIZE IS WHITE PP.

SIZE CODE.....:	-0030		-0060		-0125		-0250	
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4		250/8	
NOM. BRIM CAP. (mL):	45		76		169		318	
NECK FINISH (nn)...:	20-415		24-415		38-430		38-430	
THREAD (T.P.I.)....:	8		8		6		6	
NOM. WEIGHT W/CLOSURE (gms)....:	14.5		20		49		60	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....:	0.54 ±0.02	13.7 ±0.5	0.70 ±0.02	17.8 ±0.5	1.100 ±0.025	27.94 ±0.64	1.100 ±0.025	27.94 ±0.64
B (HT. W/CLOSURE)..:	2.50 ±0.05	63.5 ±1.3	3.25 ±0.06	82.6 ±1.5	4.29 ±0.08	109.2 ±2.0	5.72 ±0.08	145.3 ±2.0
C (HT. W/O CLOSURE):	2.39 ±0.05	60.7 ±1.3	3.15 ±0.06	80.0 ±1.5	4.12 ±0.08	104.6 ±2.0	5.55 ±0.08	141.0 ±2.0
D (SQUARE).....:	1.50 ±0.05	38.1 ±1.3	1.58 ±0.06	40.1 ±1.5	2.05 ±0.06	52.1 ±1.5	2.30 ±0.06	58.4 ±1.5
W.T.(NOM. WALL THK):	0.04	1.0	0.04	1.0	0.06	1.5	0.06	1.5

SIZE CODE.....:	-0500		-1000		-2000	
SIZE (NOM VOL mL/oz)	500/16		1000/32		2000/64	
NOM. BRIM CAP. (mL):	595		1210		2380	
NECK FINISH (nn)...:	38-430		38-430		53B*	
THREAD (T.P.I.)....:	6		6		4	
NOM. WEIGHT W/CLOSURE (gms)....:	108		170		325	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....:	1.100 ±0.025	27.94 ±0.64	1.100 ±0.025	27.94 ±0.64	1.53 ±0.03	38.9 ±0.8
B (HT. W/CLOSURE)..:	6.98 ±0.09	177.3 ±2.3	8.63 ±0.10	219.2 ±2.5	10.65 ±0.14	270.5 ±3.6
C (HT. W/O CLOSURE):	6.80 ±0.09	172.7 ±2.3	8.43 ±0.10	214.1 ±2.5	10.42 ±0.14	264.7 ±3.6
D (SQUARE).....:	2.90 ±0.08	73.7 ±2.0	3.65 ±0.07	92.7 ±1.8	4.57 ±0.10	116.1 ±2.5
W.T.(NOM. WALL THK):	0.06	1.5	0.05	1.3	0.09	2.3



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D	PSS 00-2015-XXXX-1	REV.
TL	GRP 10/23/09		1009
CK'D DM 10/22/09			

DESCRIPTION: SQUARE, NARROW-MOUTH GRADUATED LAB QUALITY BOTTLES