
CATALOG NUMBER	POR SIZE (µm)	MATERIALS		A (HEIGHT)		ØB (DIAMETER)	
		HOUSING	MEMBRANE	INCH	MM	INCH	MM
724-2080	0.8	PP	CA	0.96	24.5	1.15	29.2
727-2020	0.2	PP	NYLON	0.96	24.5	1.15	29.2
727-2045	0.45	PP	NYLON	0.96	24.5	1.15	29.2
728-2020	0.2	PP	PFFE	0.96	24.5	1.15	29.2
728-2045	0.45	PP	PTFE	0.96	24.5	1.15	29.2
724-2020	0.2	PP	SFCA	0.96	24.5	1.15	29.2
724-2045	0.45	PP	SFCA	0.96	24.5	1.15	29.2



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETORY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	APP'D		REV.
PMP	PWC	PSS 72B-XXXX	0216
CK'D	2/22/16		

DESCRIPTION:

SYRINGE FILTERS, NON-STERILE