Technical Data Sheet

Article Number	**41408002**
Description	**Glass wool, extra-fine, packed per 30 g**

Hydrolytic resistance of glass grains according to ISO 719	HGB 3
Acid class DIN 12116	3
Alkali class DIN 52322	2
Mean fiber diameter [µm]	11
Density [g/cm^3]	2.52
Lubricant add-on [%]	Max. 0.2 basing on synthetic oils
Transformation temperature [°C]	530
Softering point [°C]	680

Glaswarenfabrik Karl Hecht
GmbH & Co KG
Stettener Str. 22-24
97647 Sondheim/Rhoen

Tel: +49 9779 8080
Fax: +49 9779 808 88
info@assistent.eu
www.assistent.eu

USt. ID: USt.IdNr./VAT-No.:
DE 132198535
Amtsgericht Schweinfurt HRA 8063
Persönlich haftende Gesellschafterin:
Hecht GmbH AG Schweinfurt HRB1160

Geschäftsführer:
Steffen Wachter
Thomas Diepholtz
Einzelprokura:
Tonja Müller