Hospidex

Safety solutions




Clinisafe RECY 60l 1PC & 2PCS

PRODUCT DATA SHEET

General data	
Item code	body : HOS-60L.PB.RECY/4 - lid 1PC : HOS-PB.LID.N - lid 2PCS : HOS-PB.LIDPH.N
GTIN code	
Product Name	Clinisafe 60 liters recycled (NL/FR/EN/DE) - lid 1PC - lid 2PCS
Description	Bin 60 liters recycled plastic (yellow) - lid 1PC (white) - lid 2PCS (white)
Brand	Clinisafe®
Business Unit	Safety Solutions
Category of safety packaging	Clinical waste bin
Indication	Segregation, protection, transport, disposal of solid sharps and non sharps biohazard waste with residual liquids

Characteristics data	
Nominal volume	60 l
Total volume	57 l
Usefull volume	49 l
Maximimun weight	20 kg
Body background aspect	Yellow
Lid color	Translucent white
Label or printing color	Yellow
Inner bag color	NC
Dominant color	Yellow
Marking languages	ENG/FR/NL/D
Dimensions (assembled) in mm :	
Long	400 mm
Wide	300 mm
Height	670 mm
Aperture	1PC : 370 x 270 mm - 2PCS : 140 x 140 mm Ø
Fill line height	590 mm
Weight	body : 1626 g - lid 1PC : 360 g - lid 2PCS : 390 g
Material	Body recycled polypropylene - lid 1PC & 2PCS virgin polypropylene
Latex free	Yes
Carbon footprint reduction	27%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with container open
Temporary closure	Yes
Permanent closure	Yes
Handles	2 lateral
Accessories availables	Trolley (see brochure and compatibility tables)





Logistic data	
Carton quantity	50 units (lid 1PC only) - 25 units (lid 2PCS only)
Carton weight	19 Kg (lid 1PC) - 11 Kg (lide 2 PCS)
Carton dimensions (LxWxH)	NC
Carton volume	NC
Palet quantity	200 units (body only)
Palet weight	320 Kg (body only)
Palet dimension (LxWxH)	120 X 80 X 200 cm (body only)

Certifications data	
ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
BSI / Kitemark	No
AFNOR / NF Mark (NF302)	No
TRBA250	NC
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	No
ISO 14001	Yes
ISO 45001	No

Manufacture data	
Made in	Italy
By	2001 srl
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group