



SUPPORT AND DRAIN RACK FOR PIPETES "SOP-PIP"

Made AISI 304 stainless steel.
Holds up to 112 pipettes of different sizes.
Part No. 1001223



TILTED SUPPORT RACK "SOP-INC"

Suitable for pipettes, hydrometers and thermometers rack, made of methacrylate. 6 positions.
Part No. 1000015



KJELDAHL FLASK SUPPORT

Made of AISI 304 stainless steel with electrobright finish.
Suitable for flasks up to 500 ml.
Part No. 1000002



PIPETTE RACK

Made of AISI 304 stainless steel.
No. of places: 6 dosing pipettes of various sizes up to 28 cm in length. (See pages 220 - 221).
Part No. 1001215



SUPPORT RACK FOR PIPETTE AID "COM-PIP"

Made of AISI 304 stainless steel.
No. of places: 6 dosing pipettes of various sizes.
(See page 220-221).
Part No. 1001702




ROUND BOTTOM FLASK SUPPORT

Made of moulded rubber (E.P.D.M.). Temperature resistance up to 100 °C. Specially designed with stepped surface to support different flask diameters. From 100 ml up to 10 litre capacity.
Dimensions: 50 mm high x 160 mm Ø.
Part No. 1000863



TABLE TOP OR WALL MOUNTABLE DRAIN RACK

Made of PCV coated wire.
640 high x 500 mm wide.
With 79 positions, complete with drip tray.
Part No. 1000481



DRAIN RACK

Made of PCV coated wire.

32 positions.
300 high x 400 mm wide.
Part No. 1024132

48 positions.
500 high x 400 mm wide.
Part No. 1024148

88 positions.
620 high x 500 mm wide.
Part No. 1024188