Technical Sheet

BLUE-RAL 555

Product code: 361650

Specifications:

• **Appearance:**	Clear blue Solution
• **Absorption Maximum in Water:**	612 nm ≤ λ_1 ≤ 617 nm 648 nm ≤ λ_2 ≤ 651 nm
• **pH:**	6.5 – 7.5

Applications:

Haematology: used in Kit RAL 555, fast variation of May-Grünwald Giemsa staining, for differential staining of medullary and blood smears.

Parasitology and Mycology : used in Kit RAL 555 to detect parasites and fungi.

Histo-Cytology: used in Kit RAL 555 for fast polychromic staining of histological sections and for the detection of Helicobacter pylori on histological sections.

Updated: 07/02/2022