Specification

PBS tablets pH 7.4 (for 100 ml) A9162

Product Code:	A9162
Product Name:	PBS tablets pH 7.4 (for 100 ml)
Specifications:	: • Ready to use in minutes : • Autoclaveable pH (25°C): 7.4 ± 0.05 **Composition:** KCL: 2.7 mM NaCl: 140 mM Phosphate: 10 mM
WGK:	1
Storage:	RT
Origin:	Tablets for preparation of 100 ml phosphate buffered saline (PBS)
CS:	38221900

AppliChem GmbH
Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • **info.de@itwreagents.com** • **www.itwreagents.com**
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340