Kartell SpA DIVISIONE LABWARE	TECHNICAL DATA SHEET ITEM 996	Mod. AQ 001 Rev. 0 del 25/09/2003

TECHNICAL SPECIFICATIONS

Trade name	Pipette Trays for drawer
Use	Very useful work station holder not only for pipette.
Material	PVC Polyvinychloride

DIMENSIONAL SPECIFICATIONS

Dimension	mm 216x283
Height	mm 40
Length	mm 129,32
Overall pipettes stem diameter accomodated	7 pipettes not exceeding 20mm longitudinally
Overall pipettes stem diameter accomodated	16 pipettes not exceeding 10mm laterally

GENERAL FEATURES

CE mark	No
Autoclavability	No
Sterilization	Gas Sterilizable (Ethylene Oxyde) Chemical Sterilizable (in formalin) Microwaves Yes.
Food contact suitable	No

CHEMICAL RESISTANCE

Categories of substances	Reference	Concen. %	T. 20°	T. 40°	T.° 60
Inorganic acid	Sulphuric acid	98	C	C	C
Organic acid	Benzoic acid	100	A	B	B
Alcohol	Ethanol	100	A	B	B
Aldehydes	Acetaldehyde	100	C	C	C
Inorganic alkali	Sodium hydroxide	50	A	A	A
Organic alkali	Aniline	100	C	C	C
Ketone	Acetone	100	C	C	C
Esther	Ethyl acetate	100	C	C	C
Hydrocarbons, halogenated	Dichloroethylene	100	-	-	-
Hydrocarbons, aromatic	Benzene	100	C	C	C
Hydrocarbons, linear	Hexane	100	B	C	C
Supersolvent	Tetrahydrofuran	100	C	C	C

Legenda (not applicable if neglected)**:**

- **A: fair resistence;** exposure (30 days) to the chemical does not cause any damage.
- **B: sufficient resistence;** exposure causes damage of poor importance, which sometimes is only temponary.
- **C: poor resistence;** exposure to chemical is not allowed, or cases immediate permanent damage.

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing..

PACKAGING

Bag	Pcs. 6

Quality Assurance — Massimiliano Capitanio

Technical Direction — Guido Borona

This paper is an electronic print-out: signature is not requested