HIRSCHMANN®

Data sheet

Article: Graduated pipette with piston
Code no.: 119 01 60

Material:
AR-Glas (Soda-lime glass)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 3 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 3 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Standard:
ISO 835 (graduated pipette)

Hirschmann Laborgeräte
GmbH & Co.KG

.....
Roland Benz
Quality Manager

Eberstadt, 10.02.2016