
Technical Drawing Braided Reinforced Silicone Tubing (TBR)

Liveo Q7- 4765

Materialnumber	OD (B)	ID (A)	±	Wall (C)	±	Length
00 45 63 60 20- (60-664-20)	9,27 mm	3,18 mm	0,25 mm	3,05 mm	0,25 mm	7,5 m
00 45 63 60 30- (60-664-30)	11,35 mm	4,75 mm	0,25 mm	3,30 mm	0,25 mm	7,5 m
00 45 63 60 40- (60-664-40)	13,21 mm	6,35 mm	0,25 mm	3,43 mm	0,25 mm	7,5 m
00 45 63 60 45- (60-664-45)	15,04 mm	7,92 mm	0,38 mm	3,56 mm	0,38 mm	7,5 m
00 45 63 60 50- (60-664-50)	16,64 mm	9,53 mm	0,38 mm	3,56 mm	0,38 mm	7,5 m
00 45 63 60 60- (60-664-60)	20,32 mm	12,70 mm	0,38 mm	3,81 mm	0,38 mm	7,5 m
00 45 63 60 70- (60-664-70)	24,51 mm	15,88 mm	0,51 mm	4,32 mm	0,51 mm	7,5 m
00 45 63 60 80- (60-664-80)	27,94 mm	19,05 mm	0,51 mm	4,45 mm	0,51 mm	7,5 m
00 45 63 60 90- (60-664-90)	31,37 mm	22,23 mm	0,51 mm	4,83 mm	0,51 mm	7,5 m
00 45 63 61 00- (60-664-100)	34,54 mm	25,40 mm	0,51 mm	4,57 mm	0,51 mm	7,5 m



17.11.2021

Quality Assurance Freudenberg Medical Europe GmbH, Van Nguyen