Technical Data Sheet



VWR® 1mL Serological Pipettes

Description

1mL in 1/100mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø4.5	±0.1
B	ø2.45	±0.35
C	268.5	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: yellow
- Graduation (mL): in 1/100 mL
- Negative graduations up to -0.3 mL (total volume = 1.3 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 1mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5504	/	1mL in 1/100mL, Individually Wrapped, Paper/Plastic	Yes	1	/	500
612-5505	75816-112	1mL in 1/100mL, Bulk packed	Yes	/	25	1000
612-5506	75816-116	1mL in 1/100mL, Bulk packed	No	/	2000	2000
612-5829	/	1mL in 1/100mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	500
612-5890	76320-448	1mL in 1/100mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	100	500
/	75816-102	1mL in 1/100mL, Individually Wrapped, Paper/Plastic	Yes	1	100	500
/	76184-872	1mL in 1/100mL, Individually Wrapped, Plastic/Plastic	Yes	1	100	500

Technical Data Sheet



VWR® 2mL Serological Pipettes

Description

2mL in 1/100mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø6.3	±0.1
B	ø2.93	±0.35
C	272.0	±2.0

STERILE | R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: green
- Graduation (mL): in 1/100 mL
- Negative graduations up to -0.3 mL (total volume = 2.3 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 2mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5507	/	2mL in 1/100mL, Individually Wrapped, Paper/Plastic	Yes	1	/	500
612-5521	75816-114	2mL in 1/100mL, Bulk packed	Yes	/	25	1000
612-5522	75816-110	2mL in 1/100mL, Bulk packed	No	/	1000	1000
612-5830	/	2mL in 1/100mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	500
612-5891	76320-274	2mL in 1/100mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	100	500
/	75816-104	2mL in 1/100mL, Individually Wrapped, Paper/Plastic	Yes	1	100	500
/	76184-742	2mL in 1/100mL, Individually Wrapped, Plastic/Plastic	Yes	1	100	500

Technical Data Sheet



VWR® 5mL Serological Pipettes

Description

5mL in 1/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø8.2	±0.1
B	ø2.5	±0.35
C	341.0	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: blue
- Graduation (mL): in 1/10 mL
- Negative graduations up to -3 mL (total volume = 8 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 5mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5826	/	5mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	400
/	76320-288	5mL in 1/10mL, Bulk packed	Yes	/	25	500
/	76320-290	5mL in 1/10mL, Bulk packed	No	/	500	500
/	76320-276	5mL in 1/10mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	50	400
/	76320-286	5mL in 1/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	100	400
/	76201-710	5mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	100	400

Technical Data Sheet



VWR® 10mL Serological Pipettes

Description

10mL in 1/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø10.96	±0.2
B	ø7.95	±0.25
C	ø2.95	±0.25
D	303.4	±2.0

STERILE | R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: orange
- Graduation (mL): in 1/10 mL
- Negative graduations up to -3 mL (total volume = 13 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation $SAL10^{-6}$ (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 10mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5541	/	10mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	400
612-5542	75816-098	10mL in 1/10mL, Bulk packed	Yes	/	25	400
612-5543	75816-096	10mL in 1/10mL, Bulk packed	No	/	400	400
612-5832	/	10mL in 1/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	400
612-5893	76320-278	10mL in 1/10mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	50	400
/	75816-100	10mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	100	400
/	76184-746	10mL in 1/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	100	400

Technical Data Sheet



VWR® 10mL Serological Pipettes

Description

10mL in 1/10mL, Wide Tip


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø10.0	±0.1
B	ø7.95	±0.25
C	ø5.16	±0.25
D	346.3	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Color coded: orange
- Graduation (mL): in 1/10 mL
- Negative graduations up to -3 mL (total volume = 13 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation $SAL10^{-6}$ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 10mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5836	76184-754	10mL in 1/10mL, Bulk packed	Yes	/	50	400

Technical Data Sheet



VWR® 10mL Serological Pipettes

Description

10mL in 1/10mL, Ultrasonic Welding


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø10.0	±0.1
B	ø7.95	±0.25
C	ø3.5	±0.25
D	346.3	±2.0

STERILE | R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: orange
- Graduation (mL): in 1/10 mL
- Negative graduations up to -3 mL (total volume = 13 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 10mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
/	76501-552	10mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	100	400
/	76501-554	10mL in 1/10mL, Bulk packed	Yes	/	25	400
/	76501-556	10mL in 1/10mL, Bulk packed	No	/	400	400
/	76501-558	10mL in 1/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	100	400
/	76501-560	10mL in 1/10mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	50	400

Technical Data Sheet



VWR® 25mL Serological Pipettes

Description

25mL in 2/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø15.1	±0.1
B	ø7.95	±0.25
C	ø5.0	±0.25
D	308.5	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: red
- Graduation (mL): in 2/10 mL
- Negative graduations up to -8 mL (total volume = 33 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 25mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5544	/	25mL in 2/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	150
612-5545	75816-092	25mL in 2/10mL, Bulk packed	Yes	/	10	150
612-5833	/	25mL in 2/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	150
612-5894	76320-280	25mL in 2/10mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	50	150
/	75816-090	25mL in 2/10mL, Individually Wrapped, Paper/Plastic	Yes	1	50	150
/	76184-748	25mL in 2/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	50	150

Technical Data Sheet



VWR® 50mL Serological Pipettes

Description

50mL in 5/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø18.4	±0.1
B	ø8.0	±0.25
C	ø5.8	±0.25
D	346.6	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: purple
- Graduation (mL): in 5/10 mL
- Negative graduations up to -10 mL (total volume = 60 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 50mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5546	/	50mL in 5/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	100
612-5837	76184-756	50mL in 5/10mL, Bulk packed	Yes	/	10	100
612-5834	/	50mL in 5/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	100
612-5895	76320-282	50mL in 5/10mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	25	100
/	75816-088	50mL in 5/10mL, Individually Wrapped, Paper/Plastic	Yes	1	50	100
/	76184-750	50mL in 5/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	50	100

Technical Data Sheet



VWR® 100mL Serological Pipettes

Description

100mL in 1mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø24.3	±0.1
B	ø7.95	±0.25
C	ø5.9	±0.25
D	346.8	±2.0

STERILE | R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: pink
- Graduation (mL): in 1 mL
- Total volume = 100 mL
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 100mL Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5547	/	100mL in 1mL, Individually Wrapped, Paper/Plastic	Yes	1	/	50
612-5548	75816-086	100mL in 1mL, Bulk packed	Yes	/	10	60
612-5835	/	100mL in 1mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	50
612-5896	76320-284	100mL in 1mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	25	50
/	75816-084	100mL in 1mL, Individually Wrapped, Paper/Plastic	Yes	1	10	50
/	76184-752	100mL in 1mL, Individually Wrapped, Plastic/Plastic	Yes	1	10	50

Technical Data Sheet



VWR® 5mL Mini Serological Pipettes

Description

5mL in 1/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø9.7	±0.1
B	ø7.95	±0.25
C	ø3.5	±0.25
D	234.95	±2.0

STERILE | R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: blue
- Graduation (mL): in 1/10 mL
- Negative graduations up to -3 mL (total volume = 8 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 5mL Mini Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5523	/	5mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	400
612-5539	75816-106	5mL in 1/10mL, Bulk packed	Yes	/	25	500
612-5540	75816-108	5mL in 1/10mL, Bulk packed	No	/	500	500
612-5831	/	5mL in 1/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	/	400
612-5892	/	5mL in 1/10mL, Individually Wrapped, Plastic/Plastic, Triple bag	Yes	1	50	400
/	75816-094	5mL in 1/10mL, Individually Wrapped, Paper/Plastic	Yes	1	50	400
/	76184-744	5mL in 1/10mL, Individually Wrapped, Plastic/Plastic	Yes	1	50	400

Technical Data Sheet



VWR® 10mL Mini Serological Pipettes

Description

10mL in 2/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø15.1	±0.1
B	ø7.95	±0.25
C	ø5.0	±0.25
D	234.95	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: orange
- Graduation (mL): in 2/10 mL
- Negative graduations up to -10 mL (total volume = 20 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 10mL Mini Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5827	/	10mL in 2/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	150
/	76204-736	10mL in 2/10mL, Individually Wrapped, Paper/Plastic	Yes	1	50	150
612-7338	/	10mL in 2/10mL, Bulk packed	Yes	/	25	150

Technical Data Sheet



VWR® 25mL Mini Serological Pipettes

Description

25mL in 5/10mL


Purpose

Use for liquid transfer

Materials

Pipet: PS (Polystyrene), Compliance with USP Class VI — Color: Clear

Filter: PP (Polypropylene) — Color: White



Unit: mm

	Dimensions	Tolerance
A	ø18.4	±0.1
B	ø8.0	±0.25
C	ø5.8	±0.25
D	234.95	±2.0

STERILE R

Features

- Temperature: -20°C to +50°C
- Available in individually packed in peel-to-open wrap
- Color coded: red
- Graduation (mL): in 5/10 mL
- Negative graduations up to -10 mL (total volume = 35 mL)
- Graduations are calibrated for accurate dispensing to within ±2%
- Autoclavability: No
- DNase/RNase free, Non-pyrogenic, Non-cytotoxic
- Latex Statement: this product is latex free
- Meet the requirement of BSE/TSE
- Sterilized by irradiation SAL10^{-6} (ISO11137)
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 25mL Mini Serological Pipettes

VWR EU Cat. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/PK	Qty/Case
612-5828	/	25mL in 5/10mL, Individually Wrapped, Paper/Plastic	Yes	1	/	100
/	76204-738	25mL in 5/10mL, Individually Wrapped, Paper/Plastic	Yes	1	25	100
612-7337	/	25mL in 5/10mL, Bulk packed	Yes	/	25	100