PRODUCT DATA SHEET

MINI-SPIKE V (4550560)



Attribute Name	Attribute Value
Brand name	Mini-Spike®
Product name (HC)	Mini-Spike® V
Product type	Dispensing pin
Application	Withdrawal and Injection
Type of filter	Air filter
Position of filter	Injection port
Pore size (Air filter)	0.45 µm
Color	green
Type of packaging	Blister
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com

Attribute Name	Attribute Value
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Disposable / Reusable	disposable
Sterility	sterile
Type of Sterilization	Gamma