

Product Specifications

Product **25862**
Buffer solution pH 5,Acetate buffer

General Product Data

Version	2
CAS No	64-19-7
Molecular weight	0.00
Molecular formula	
Linear formula	

Product Specifications

Appearance (Color)	Clear colorless
Appearance (Form)	Solution
pH	4.98 to 5.02 (at 20°C) composition: 0.7 ml glacial acetic acid and 0.6 g potassium hydroxide per liter

Product Specifications are subject to amendment and may change over time.

Acros Organics
ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.acros.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329