Certificate

Filter paper MN 625

REF: 543024

Product description

Slow filtration, for fine precipitates

typical data	unit	specifies	based on
basis weight	g/m²	130	DIN EN ISO 536
thickness	mm	0.26	DIN EN ISO 12625-3
filtration speed	sec	30	DIN 53137
surface		smooth	

Confirmation

Hereby we confirm that the above-mentioned product has successfully passed our quality control system in accordance with EN ISO 9001 and meets the specific quality criteria.

This document has been produced electronically and is valid without a signature.

Management System EN ISO 13485:2016 ISO 9001:2015
TÜVRheinland CERTIFIED
www.tuv.com ID 9000008338

MACHEREY-NAGEL GmbH & Co. KG
Valencienner Str. 11
52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com