Hospidex

Datasheet

SHA-51125610



5th GEN. 7 liters wide with Dual Protect Pro

General data

Item code	SHA-51125610
GTIN code	5056271508948
Product Name	Sharpsafe® 5th GEN. 7 liters wide with Dual Protect Pro
Description	Sharps container 7 liters wide componed by 1 body + 1 lid to assemble
Brand	Sharpsafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste

Manufacture data

Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group

Hospidex nv

Grijpenlaan 23, B-3300 Tienen • info@hospidex.eu • T. +32 (0)16 77 89 31

www.hospidex.eu

Characteristics data

Nominal volume	7 l
Total volume	6,96 l
Usefull volume	5,85 l
Maximum weight	1,40 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
	Dimensions (assembled) in mm
Long	205 mm
Wide	170 mm
Height	342 mm
Opening	65 x 61 mm
Fill line height	287 mm
Weight	405 g
Material	Polypropylene recycled 100% (body) + Polypropylene virgin 100% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction footprint Vs virgin material	67%
Autoclavable	Yes
Autoclavability conditions	134° during 18 mn with lid open
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availibility Accessories	Stands, trolley brackets (see brochure and compatibility tables)

Logistic data

Quantity per carton	48
Carton weight	20,4 Kg
Carton dimensions (LxWxH)	595 x 356 x 712 mm
Carton volume	0,151m3
Pallet quantity	480 containers
Pallet weight	204 Kg
Pallet dimension (LxWxH)	1200 x 800 x 2298 cm

Certifications data

ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark - (NF 302 in progress)
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
	Quality management system
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes

Hospidex nv
Grijpenlaan 23, B-3300 Tienen • info@hospidex.eu • T. +32 (0)16 77 89 31
www.hospidex.eu

SS-DS-H0261-SharpsafeSGEN7Large-EN 202509