DURAN
WHEATON
KIMBLE

Datenblatt - *Data Sheet*

Article:

Tube, Centrifuge, OIL, WEATHERING, 100 mL

Can be used in the determination of residues in Liquified Petroleum (LP) gases, ASTM Method D2158. Calibrated to contain. Different graduations than 45240. Scale and legend are durable white ceramic enamel. Top is tooled to accept a Snap Cap (28150R-6).

Article No:

45241-100


Material:

Glass: **Boro 3.3 – ASTM E438**

Autoclavable

Dimension:

Height H [mm]	max. 203
Diameter D [mm]	Ø37
Diameter Di [mm]	Ø17
Diameter Da [mm]	Ø28
Volumina Tolerances [mL]	0 to 0.1 - ±0.02
	Above 0.1 to 0.30 - 0.03
	Above 0.3 to 1 - 0.05
	Above 1 to 3 - 0.10
	Above 3 to 5 - 0.20
	Above 5 to 10 - 0.50
	Above 10 - 1.00



Subject to technical modifications!