TECHNICAL DATA SHEET

29076983

CONT-DAIRY-SNAPSEAL--PP---CLEA--8,0--

8 OZ SNAP SEAL VIAL

Single Use Vial for DAIRY,
Food Sampling

Snap Seal Closure



Applications

Intended Fields of Use

- Field collection of farm milk and related dairy products
- Quality control sample storage containers and food packaging vials

Options

- Can be supplied with label (barcode…)
- Standard vial is colorless – colors upon request

Characteristics

Characteristics	Typical Values
Internal Volume	240ml
Filling Level Line	160 ml
Vial Overall Height	76 mm
Top Body Diameter	77 mm
Bottom Body Diameter	65 mm
Vial Weight	30,6 grams

Delivery

Qty Per Box	120 units
Qty Per Pallet	3600 units

Benefits

Vial Design

- User-friendly hinged vials concept
- Patented “in mold closing” process
- Translucent in color for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use
- Snap Seal closure secures samples
- Allows the use of gloves in cold environment

Regulatory and Quality

- ISO-9001 certified company
- Complies with Regulation (EC) No 1935/2004

Please contact us for official regulatory submission documents.

(Rev. 06/2024)