multiroir
LA MEILLEURE FAÇON DE RANGER

SKU : MULC50206BAC400ET

BOX STORAGE DIVIDABLE - SAN 400E

Model : BOX STORAGE DIVIDABLE



The SAN 400 E divisible bin is an ideal storage solution for laboratory technicians, meeting the specific demands of laboratory environments. Its complete transparency allows for immediate identification of samples and reagents, optimizing the efficiency of analytical procedures and reducing the risk of errors. Designed to seamlessly integrate into a laboratory workspace, this bin can be safely placed in storage cabinets, on workbenches, or shelves, thanks to its balanced and stable design. The recessed handle makes handling easy, even with laboratory gloves, ensuring the safe transport of sensitive samples or chemicals.

The optimized dimensions of 400 x 91 x 81 mm make the SAN 400 E bin particularly suitable for storing a variety of laboratory equipment, from test tubes to pipettes and Petri dishes. Its ribbed bottom is a major advantage for laboratory technicians, facilitating the retrieval of small items such as microscope slides or tube caps, thus saving time and improving work precision. Additionally, the secure stackability of these bins allows for optimal organization of laboratory space, a crucial aspect in often constrained environments. Made of SAN, this bin offers excellent resistance to chemicals commonly used in laboratories, ensuring its durability and integrity under intensive use conditions.

PRODUCT ADVANTAGES

- **Monobloc construction :** Made from a single piece, this design ensures unparalleled reliability.
- **Stackable:** Its design allows the bin to be stably stacked with other models.
- **Robust :** Entirely made of SAN, the bin offers excellent resistance to chemical agents, oils, and temperature variations ranging from -20°C to +80°C.
- **Easy handling :** The recessed handle ensures easy manipulation of the bin, even when loaded.

SPECIFICATIONS

SKU	MULC50206BAC400ET
External dimensions	400 x 91 x H81 mm
Internal dimensions	365 x 80 x H70 mm
Wall/Base type	Solid / Ribbed
Handle Type	Recessed
Volume	2 L
Material	SAN
Color	Transparent
Weight	250 g
Stackable	Yes
Operating temp.	-20°C / +80°C
Shipping time	5 working days

Multiroir 10 route de Brie Comte Robert 94520 Périgny-Sur-Yerres FRANCE — **Téléphone** +33 (0)1 45 10 15 20 — **Fax** +33 (0)1 45 10 15 95
info@multiroir.com — www.multiroir.com