

Model:
PLED10 - 631-1086

Typology:
ILLUMINATOR

Description:
Double flexible arm (L = 60 cm) with double P-LED³

Chromaticity coordinates (x,y):	0.320, 0.320 (pure white)
Electrical power	3,6W
Optical power:	Luminous flux: 400 lumen Illuminance (at 10 cm distance): 170000 lux
Aperture angle	25°
Power supply:	Brightness control potentiometer. External power supply provided. Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA