
Codes: 300250, 300259, 300261, 300252, 300251, 300253, 300254, 300263, 300250.1
Sterile swab in round tube

Mod.	Code	Description	Case Quantity	Case Weight	Case Volume	Cases per pallet
B	300250	Wood + cotton	4 x 500	15.00	0.065	24
B	300259	Wood + viscose	4 x 500	14.00	0.070	24
A	300261	Snappable PS + cotton	4 x 500	14.20	0.070	24
A	300252	Snappable PS + viscose	4 x 500	14.20	0.070	24
C	300251	Aluminium + cotton	4 x 500	14.00	0.070	24
C	300253	Aluminium + viscose	4 x 500	13.90	0.070	24
D	300254	Twisted aluminium + viscose	4 x 500	13.90	0.070	24
A	300263	Snappable PS + polyester	4 x 500	14.20	0.070	24
B	300250.1	Wood + pure cotton	4 x 500	14.00	0.070	24

Directive 93/42/EEC. Medical devices



Elaborated	Verified	Authorised
Gemma Riba Marketing Dept.	Anna Mir Technical, Development and Quality Director	Mònica Torras Sales & Marketing Dir.

SPECIFIC USE AND DESCRIPTION

Sterile swabs for clinical use, designed for the collection of biological samples and for their subsequent transport to the laboratory for microbiological analyse.

Mainly used for the preparation of smears (stained and then observed through a microscope). Also useful for isolations in culture media.

Another remarkable use is for surface spreading techniques as required in antibiotic susceptibility testing for instance.

Another remarkable use is based on its application to plant plates by the technique of spreading on the surface, as in the case of susceptibility.

It is also useful for multiple uses as cleaning biological cavities or cleaning inaccessible spots, among others.

Expiry date: 48 months from sterilisation date (sterilisation by ethylene oxide). Certificates of sterilisation are available by request.
The swabs are considered as surgically invasive products for a transient use. And therefore, swabs are ISO 13485 manufactured and have the . Copies of both certificates are available by request.

The CE MDD, Class IIa Mark ensures that the facilities and the dynamic work of the company are subjected to rigorous internal and external quality controls. In our case, the evaluation of infrastructure and external monitoring is made by the AEMPS CE 0318 (Spanish Agency of Medicines and Medical Devices).

COMPONENTS

- A 13 x 165 mm sterile round bottom tube with cap, made of non-distorting polypropylene and suitable for food use.

- A swab consisting of a shaft (made of wood, aluminium or snappable PS) equipped with a safety cap that perfectly seals the tube (cap made of red polyethylene).

- A label that seals the cap of the tube. The label provides space for patient's name, date and hour of sample collection, number, Doctor, specimen nature and hospital name. The expiry date, batch number, product description and brand are also indicated.

Elaborated	Verified	Authorised
Gemma Riba Marketing Dept.	Anna Mir Technical, Development and Quality Director	Mònica Torras Sales & Marketing Dir.

DIMENSIONS

Code	Shaft + tip	Shaft length mm	Shaft Ø mm	Tip Ø mm	Tube int. Ø mm	Tube ext. Ø mm	Tube + cap length mm
300250	Wood + cotton	150	2.5	5	10.94	12,85	145 + 20
300259	Wood + viscose	150	2.5	5	10.94	12,85	145 + 20
300261	Snappable PS + cotton	150	2.5	5	10.94	12,85	145 + 20
300252	Snappable PS + viscose	150	2.5	4	10.94	12,85	145 + 20
300251	Aluminium + cotton	150	0.9	1.5	10.94	12,85	145 + 20
300253	Aluminium + viscose	150	0.9	1.5	10.94	12,85	145 + 20
300254	Twisted alumin. + viscose	147	0.6	3	10.94	12,85	145 + 20
300263	Snappable PS + polyester	150	2.5	5	10.94	12,85	145 + 20
300250.1	Wood + pure cotton	150	2.5	5	10.94	12,85	145 + 20

Elaborated	Verified	Authorised
Gemma Riba Marketing Dept.	Anna Mir Technical, Development and Quality Director	Mònica Torras Sales & Marketing Dir.