

For prices and availability, please log in or wait for the validation of your account

Log in to add product to cart

5150T060-qp

Range : -10+60
Unit : °C
Div. : 0.5
Accuracy : ± 0.8
Immersion : TOT
Length mm : 305
Package : QP

They are manufactured using yellow Sodocalcique glass with a diameter 6-7 mm.
Mercury free : red liquid filling.
They feature a special ring/hook that can be clipped on the glass, patented by Alla France.

They possess an expansion chamber to support temperatures higher than those of the thermometer.

The graduation is printed and fused in the glass with a special process. The graduation resists to all chemical products except fluoridric acid and hot soda (to a lesser extent) because these 2 products attack the glass. Even if you can feel a relief when touching the graduation, the thermometer is not engraved.

Available in partial or total immersion.
For partial immersion (76mm usually), there is a line printed around the thermometer at 76mm with a clear mention: "Imm. 76mm" enables to immerse the thermometer to the correct depth..
Individually packed in recyclable sturdy plastic concertina box s, ideal for storing the instrument when not in use.