

UV LAMPS

Lamps 6W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-6.LC	365 / 254	145 x 48	2 x 6	UV Fluorescence
VL-6.L	365	145 x 48	1 x 6	UV Fluorescence
VL-6.C	254	145 x 48	1 x 6	UV Fluorescence

VL-206.BL	365	No Filter	2 x 6	Black Light
VL-206.G	254	No Filter	2 x 6	Germicidal

VILBER
info@vilber.com - www.vilber.com



UV LAMPS

Lamps 8W

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-8.LM	365 / 312	220 x 48	2 x 8	UV Fluorescence
VL-8.LC	365 / 254	220 x 48	2 x 8	UV Fluorescence
VL-8.MC	312 / 254	220 x 48	2 x 8	UV Fluorescence
VL-8.L	365	220 x 48	1 x 8	UV Fluorescence
VL-8.M	312	220 x 48	1 x 8	UV Fluorescence
VL-8.C	254	220 x 48	1 x 8	UV Fluorescence
VL-208.BL	365	No Filter	2 x 8	Black Light
VL-208.G	254	No Filter	2 x 8	Germicidal

VILBER
info@vilber.com - www.vilber.com



Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-215.L	365	295 x 76	2 x 15	UV Fluorescence
VL-215-M	312	295 x 76	2 x 15	UV Fluorescence
VL-215-C	254	295 x 76	2 x 15	UV Fluorescence
VL-115.L	365	295 x 66	1 x 15	UV Fluorescence
VL-115-M	312	295 x 66	1 x 15	UV Fluorescence
VL-115-C	254	295 x 66	1 x 15	UV Fluorescence
VL-315-BL	365	No Filter	3 x 15	Black Light
VL-315-G	254	No Filter	3 x 15	Germicidal
VL-215-BL	365	No Filter	2 x 15	Black Light
VL-215-G	254	No Filter	2 x 15	Germicidal
VL-115-BL	365	No Filter	1 x 15	Black Light
VL-115-G	254	No Filter	1 x 15	Germicidal



Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-320.BL	365	No Filter	3 x 20	Black Light
VL-320.G	254	No Filter	3 x 20	Germicidal
VL-220.BL	365	No Filter	2 x 20	Black Light
VL-220.G	254	No Filter	2 x 20	Germicidal
VL-120.BL	365	No Filter	1 x 20	Black Light
VL-120.G	254	No Filter	1 x 20	Germicidal



Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-320.BL	365	No Filter	3 x 30	Black Light
VL-320.G	254	No Filter	3 x 30	Germicidal
VL-220.BL	365	No Filter	2 x 30	Black Light
VL-220.G	254	No Filter	2 x 30	Germicidal
VL-120.BL	365	No Filter	1 x 30	Black Light
VL-120.G	254	No Filter	1 x 30	Germicidal



Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-320.BL	365	No Filter	3 x 40	Black Light
VL-320.G	254	No Filter	3 x 40	Germicidal
VL-220.BL	365	No Filter	2 x 40	Black Light
VL-220.G	254	No Filter	2 x 40	Germicidal
VL-120.BL	365	No Filter	1 x 40	Black Light
VL-120.G	254	No Filter	1 x 40	Germicidal



UV LAMPS

Irradiation Surface For A Height Of About 2m

Power (W)	1 Tube	2 Tubes	3 Tubes
15	1.5 x 2m	1.5 x 3m	1.5 x 4m
20	2m x 2m	2m x 3m	2m x 4m
30	2.5m x 2.5m	2.5m x 3.5m	2.5m x 4.5m
40	3m x 3m	3m x 4m	3m x 5m

VILBER
info@vilber.com - www.vilber.com