KIMTECH*

Product Data Sheet for Kimtech® Pure W4 Wipers - Individual / White - 7646

Product Information

Product Long Code	07646021
Customer Facing Code	7646
Description	Kimtech® Pure W4 Wipers - Individual / White
PDS Reference	24/01/2023

Packaging Configuration

Case Contents	7646: 5 Packs x 100 Sheets = 500 Sheets
Outer Pack Material	Corrugated

Icon	Unit	Dimension					Weight(kg)
	Pack	23.00	x	23.00	x	6.60	Length x Width x Height(cm)	0.38
	Sheet	22.80	x	22.80			Length x Width(cm)	

Product Description

Kimtech® Pure W4 Wipers are a range of high-performance wiper solutions, suitable for ISO Class 4 or higher clean-room environments. They are perfect for dry wiping of critical equipment and surfaces in pharmaceutical, high-tech biotechnology, medical device and food processing industries. These 100% Polypropylene white base sheets are acid, base and solvent resistant. They are also thermally bonded with no glue or binding and have straight, clean, knife-cut edges to minimise lint. W4 Wipers are doubled-bagged with a case liner. Individual sheet size 22.8cm (L) x 22.8cm (W). 5 packs x 100 cloths, white (500 total)

Product Features

- 5 packs x 100 white cloths (500 total)
- Low linting, 100% Polypropylene white base sheets
- Perfect for pharmaceutical, high-tech biotechnology, medical device and food processing industries
- Suitable for ISO Class 4, or higher, clean-room environments
- Acid, base and solvent resistant
- Doubled-bagged with a case liner

Product Specifications

Fold Type	Flat	**Impregnated**	No
Bar Code (Case)	10036000333906	**Bar Code (Pack)**	036000333909
Virgin Fibre(%)	100	**Technology**	Meltblown
Colour	White		

Product Data Sheet for Kimtech® Pure W4 Wipers - Individual / White - 7646

International Standards and Symbols

Solution

Safer Workplace

We know staff safety is your top priority. You can rely on our solutions to protect them, their environment and your important work.



Productive Workplace

The success of your business hinges on your work efficiency. That's why our solutions are designed to help improve your productivity.



Trade

ISO 14001

ISO 14001 is an international environmental standard. Designed to help organisations ensure that operations comply with environmental laws, and; major environmental risks and liabilities are properly identified, minimised and managed.

ISO 9001

This gives the requirements for quality management systems, is now firmly established as the globally implemented standard for providing assurance about the ability to satisfy quality requirements and to enhance customer satisfaction in supplier-customer relationships.

Packaging

Unit Of Sale	Case	**Net Weight(kg)**	2.32
Max Weight(kg)	2.78	**Case Dimension (LXWXH)(cm)**	23.8 x 23.8 x 34.9
Inner Pack Material	Mixed Plastic		

Pallet Type	LxWxH (m)	Cases per Layer	Layers per Pallet	Cases per Pallet
S2	1.20 x 1.00 x 1.21	20	3	60

Information

It is the employer's responsibility to assess the risk of the task to be undertaken and determine the correct choice of wipers for the task. The manufacturer, Kimberly-Clark, does not accept any responsibility for the incorrect choice or misuse of wipers shown in this brochure. All care has been taken to ensure that the information contained herein is as accurate as possible at the time of publication, however errors may occur and legislation concerning personal protective equipment is under constant review and may change in the lifetime of this brochure. Accordingly, the specifications for the products may be subject to change. We would advise you to contact INFOFAX if you have any queries concerning the products shown or the suitability of such products for the particular task.Always dispose of used protective equipment in a safe and appropriate manner in accordance with European, National and Local environmental regulations

Websites	www.kcprofessional.co.uk
Customer Service Number	For more information please contact Customer Service during regular office hours-0870 551 044
Country of Origin	UNITED STATES