Sodium Hypochlorite Solution (5% Available Chlorine)
BAKER ANALYZED® Reagent

Material No.: 9416-01
:
Revision No.: 0

Product Specification

Test	Specification
Appearance (Clear, light yellow liquid)	Passes Test
Available Chlorine (when packaged)	5.0 - 6.0 %
Phosphate (PO_4)	≤ 0.0005 %
Calcium (Ca)	≤ 0.001 %

For Laboratory,Research,or Manufacturing Use

Packaging Site: Paris Mfg Ctr & DC

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700
Avantor Performance Materials, LLC
100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700