Technical Data Sheet

Multi-Use Box - 1.75 L





Description: Transport and storage containers, Multi-Use Box

Purpose: Made of transparent virgin polypropylene which is food grade. Multi-use boxes are strong, transparent plastic containers that can be used to package, store and transport a wide variety of products, including foodstuffs.

Materials

Polypropylene





Unit: mm

	Dimensions	Tolerance
A	210	± 2.0
B	140	± 2.0
C	100	± 2.0

Features

- Transparent colour makes the contents easily visible
- Hygienic smooth walls for easy cleaning
- Food grade, freezer and dishwasher safe
- Suitable for –25 to +95 °C storage
- Stackable when lid is in place, nestable without the lids when not in use
- 90% space saving when empty
- Lid is secured in place with the clip handles

Multi-Use Box - 1.75 L

VWR Cat. No.	Description	Outer dim. W×D×H	Internal dim. at top W×D	Internal dim. at bottom W×D	Qty/Pk
216-2378	MULTI USE BOX - 1.75L	210 x 140 x 100 mm	170 x 120 mm	165 x 110 mm	10