Product Specifications

Product **123990010**
Diethyl ether,99+%,pure,stabilized with BHT

General Product Data

Version	6
CAS No	60-29-7
Molecular weight	74.12
Molecular formula	C4 H10 O
Linear formula	(CH3CH2)2O
Flash point (°C)	-45

Product Specifications

Appearance (Form)	Clear liquid
Color scale	=<10 APHA
Infrared spectrum	Conforms
Refractive index	1.3510 to 1.3540 (20°C, 589 nm)
GC	>=99.0 %
Water	=<0.2 % (Coulometric)
Peroxides	=<0.001 %
Residue after evaporation	=<0.002 %
Stabilizer	5 to 10 ppm BHT

Product Specifications are subject to amendment and may change over time.

Acros Organics ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: https://www.thermofisher.com

1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329