44457 n-Heptane, anhydrous, 99+%, packaged under Argon in re-sealable ChemSeal™ bottles



Product Number:	44457
CAS number:	142-82-5
MDL number:	MFCD00009544
Molecular formula:	C_7H_{16}
Linear formula:	$CH_3(CH_2)_5CH_3$
Molecular weight:	100.21

Product Specification

Appearance (Color):	Colorless
Form:	Liquid
Assay (GC):	≥99%
Water Content (Karl Fischer Titration):	≤10 ppm
Residue on Evaporation:	≤5 ppm

Date of Print:	May 22, 2022
Version:	1

Product Specifications are subject to amendment and may change over time.