SAMPLING SYSTEMS	**Product Specification**	Specification Number	8045A-1000S
		Revision	06

Part No.	8045A-1000S
Description	SteriWare PharmaScoop® – 1000ml
Use:	Single use, disposable scoop
Material of Construction:	High Impact Polystyrene (HIPS) – rigid, white - FDA acceptable CFR 177.1640 - EU Regulations 10/2011 compliant



Method of Construction:	Single piece, injection moulded
Overall Length:	324mm
Handle Length:	149mm
Maximum Width:	150mm
Nominal Weight:	169g
Assembly & Packing Environment:	Class 100,000 Medical Cleanroom
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Samplers per Box:	20 (Supplied as 5 boxes of 20 pcs)
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

SteriWare and PharmaScoop are registered Trademarks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature:	Name: Jonathan Fenner Position: Director Date: 23-Jul-12

Sampling Systems Ltd.
4-5 Forge Mills Park, Station Road, Coleshill. Warwickshire. B46 1JH. UK.
Tel: +44(0)1675 466992 Fax: +44(0)1675 466994
Email: info@sampling.com Web: www.sampling.com