
Technical Product Data

Technical Data Sheet | Labcon Product Number 3813-870-008-9



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

SuperClear® 5mL Specimen Collection and Transport Tubes, Bulk

Fluid Measure	mL	**Freezable**	Yes
Fluid Volume	5	**Freeze Temp**	-80°C
Height (mm)	60.2	**Autoclavable**	121°C 15min
Width (mm)	15.2	**Sterile**	No
Length (mm)		**Sterile Processing**	NA
Pack Qty	500	**Centrifugation**	Yes
Case Qty	1000	**Centrifugation Speed (RCF)**	12000
Package Type	Bulk	**Shelf Life (Months)**	48
Material	Polypropylene	**Certification**	ISO 9001
Cap Material	Polypropylene	**Regulatory**	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

IATA 95kPa certified to conform with USDOT HMR, 49 C.F.R., Parts 171-180 and IATA DGR UN3373 PI650 mailing tube standards and packaging instructions for diagnostic and infectious specimens.

This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is **© Copyright 06-Sep-2023** . All Rights Reserved. Technical Data Sheet Version 1 Rev 1/2023