DURAN
WHEATON
KIMBLE

# Datenblatt - *Data Sheet*

| | |
|---|---|
| *Article* | **Test Tube 6x50**<br>Boro 5.1, without rim round bottom |
| *Article Number* | **73500 – 650** |
| *Material:* | **Borosilicate Glass 5.1**<br><br>•Water resistance tested to ISO 719 (at 98 °C): HGB 1<br>•Alkali resistance tested to ISO 695: class A2<br>•Acid resistance tested to ISO 12116: S1<br>•Maximum using temperature long-time load 400°C (*)<br>•Maximum using temperature Short-time load 500°C (*)<br>•Autoclavable |
| *Dimension/Standard:* | **Gemäß ISO 4142 / refering to ISO 4142** |

| Over all length H | Outside diameter D | Wall W |
|---|---|---|
| 50 mm | 6 mm | 0,5 mm |

