Nasco SAMPLING
901 Janesville Avenue, Fort Atkinson, WI 53538-0901

Whirl-Pak® Sample Bag
Material Specification Sheet
B01065WA

Size: 4 ½" W x 9" L (11 ½ x 23 cm).

Capacity: 18 oz. (532 ml).

Film: Blend of low-linear and low-density virgin polyethylene, 2.5 mil (.064 mm) thick, FDA approved for food contact.

Write-on Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals.

Closing Tab: Aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Sterilized with ethylene oxide gas, documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality system certified to **ISO9001 & ISO13485.**

Packaging: 500 bags/box in a liner bag twist-tied shut.

Condition of Use: Not to be used at temperatures above 180° F (82° C).