




Containers

Bottles, tubes, jerry cans

Narrow-necked bottles for universal use | easy to handle and squeeze, reliable outflow





Universal bottles for storing, packing, lab, taking samples, sampling and transport of liquids.

- With screw cap
- LDPE, transparent
- With wash bottle closure
  see Page 174
- With dropper cap
  see Page 176

Bulk prices

Narrow-necked bottles for universal use

Capacity ml	Container dia. mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.
10	26	47	9	14		100		0302-0010 N3
20	32	58	9	14		100		0302-0020 N3
30	35	67	9	14		100		0302-0030 N3
50	40	86	13	18		100		0302-0050 N3
100	45	106	13	18		100		0302-0100 N3
200	58	134	13	18		25		0302-0200 N3
250	60	140	18	25		50		0302-0250 N3
500	75	173	18	25		50		0302-0500 N3
1000	95	213	21	28		50		0302-1000 N3
2000	120	263	21	28		50		0302-2000 N3
3000	139	306	24	32		20		0302-3000 N3
5000	163	360	31	40		10		0302-5000 N3

Round bottles with snap closure



Handy round bottles with snap-action hinged lid, especially suited for cosmetic liquids, emulsions, creams and oils. Sampling possible per portion.

- HDPE transparent
- Snap closure with hole Ø 3 mm
- Cost-effective and shapely

Bulk prices

Round bottles with snap closure


Capacity ml	Height with cap mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
100	115	22		50		0308-0100	N3
250	143	22		50		0308-0250	N3
500	182	25		50		0308-0500	N3
1000	221	25		72		0308-1000	N3

Tall-form bottles



Universal bottle for storing, packing, lab, sampling.

- LDPE, transparent
- Leakproof
- With screw cap
- Matching dropper cap is available as accessory



Tall-form bottles

Capacity ml	Dia. mm	Height mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
100	45	98	13	18		100		0314-0100	N3
250	61	128	18	25		50		0314-0250	N3
500	77	162	18	25		50		0314-0500	N3
1000	95	198	21.6	28		50		0314-1000	N3

Accessories

Description	Price €	Item no.	
① Dropper cap fits for thread dia. 18 mm		0300-3018	N3
① Dropper cap fits for thread dia. 25 mm		0300-3025	N3
① Dropper cap fits for thread dia. 28 mm		0300-3028	N3

The purity of disposable products and ecological sustainability combined perfectly!









With our new LaboPlast® Bio and SteriPlast® Bio product lines, we offer you products that features all the properties for high-quality sampling.
The improved specifications in this area can in part only be realised with disposable products. However to take into account environmental awareness and to safeguard ecologically sustainable trade, we now also offer these disposable products made of bio-plastic.

LaboPlast® bio-products are only manufactured from renewable raw materials. SteriPlast® bio-products are also sterilised by gamma rays. To create a completely sustainable product, the individual packaging is also made of bio-PE. This is 100 % recyclable and is manufactured without fossil raw materials. This also involves a CO_2 saving as well as a reduction of the greenhouse effect.



Sugarcane
100 % renewable raw material



100 % recyclable



I'm green™

Green PE

Reduces CO_2 (2.1 t/t PE)
Reduces the greenhouse effect

LaboPlast Bio

SteriPlast Bio

Optimum product properties

- A perfect match for sampling purposes
- With EU foodstuffs and FDA approval
- Mainly manufactured, assembled and packed in Bürkle's own clean room
- SteriPlast® Bio sterilised by gamma rays is also available

Environmentally friendly and sustainable

- Product and packaging made of bio-plastic
- Made of sustainable raw materials, dispensing with fossil materials
- CO_2 saving and thus reduction of the greenhouse effect
- 100 % recyclable
- Sugarcane, 100 % renewable raw material – reduces CO_2 (2.1 t/t PE) – reduces the greenhouse effect

Narrow-necked bio-bottles PE

The narrow-necked bio-bottles are manufactured from green PE biopolymer. This material has the same outstanding properties as conventional PE and can be recycled. However it is manufactured from renewable raw materials instead of fossils. Green PE therefore helps to reduce CO_2 emissions and protect the environment.

- With screw connection
- Manufactured from renewable raw materials
- Green LDPE



LaboPlast Bio

Narrow-necked bio-bottles PE

Capacity ml	Container dia. mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.
10	26	47	9	14		100		0303-0010 N3
20	32	58	9	14		100		0303-0020 N3
50	40	86	13	18		100		0303-0050 N3
100	45	106	13	18		100		0303-0100 N3
250	60	140	18	25		100		0303-0250 N3
500	75	173	18	25		100		0303-0500 N3
1000	95	213	21	28		100		0303-1000 N3

Laboratory bottle PET sterile



Sterile laboratory bottle for use in laboratories and chemical handling areas. Good product presentation due to crystal-clear PET. Space-saving and absolutely leakproof. Ideal for storing and packing liquids.

- Sterile SAL 10^{-3} (according to ISO 11137)
- PET crystal-clear
- With tamper-proof closure
- With moulded-in raised graduation

Laboratory bottle PET sterile

Capacity ml	L x W mm	Height with cap mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.	
125	53 x 53	105	32		336		0370-0125	N3
250	64 x 64	132	32		144		0370-0250	N3
500	76 x 76	173	32		120		0370-0500	N3
1000	96 x 96	212	32		48		0370-1000	N3

Laboratory bottles with graduation PC



High-quality laboratory bottle made of polycarbonate for universal use. Shatterproof at room temperature. Easily operated screw thread, safe pouring, with moulded-in raised contents graduation.

- Polycarbonate crystal-clear, autoclavable
- With screw cap PP
- Guaranteed leakproof
- With moulded-in raised graduation

Laboratory bottles with graduation PC

Capacity ml	L x W mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
125	55 x 55	109	28	37		6		0316-0125	N2
250	61 x 61	145	28	37		6		0316-0250	N2
1000	95 x 95	218	28	37		4		0316-1000	N2
2000	116 x 116	275	39	53		6		0316-2000	N2

Narrow-necked square bottle



Cost-effective square bottle for packing and storing of liquids, powders and granules.

- HDPE, transparent, shatterproof
- Special labelling area for labels or manual lettering
- Leakproof

Narrow-necked square bottle

Capacity ml	L x W mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
20	27 x 27	51	13	18		100		0331-0020	N3
50	33 x 34	77	13	18		100		0331-0050	N3
100	43 x 43	78	17.5	25		100		0331-0100	N3

Narrow-necked square bottle, screw cap



Cost-effective square bottle for packing and storing of liquids.

- HDPE, transparent, shatterproof
- Special labelling area for labels or manual lettering
- Suitable as a disposable bottle
- Leakproof
- With screw cap

Bulk prices

Narrow-necked square bottle with screw cap

Capacity ml	L x W mm	Height with cap mm	Thread dia. mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.	
100	41 x 33	113	18	13		100		0332-0100	N3
250	51 x 43	160	25	17.5		50		0332-0250	N3
500	68 x 53	182	25	17.5		50		0332-0500	N3
1000	89 x 66	230	25	17.5		50		0332-1000	N3

Accessories

Description	Price €	Item no.	
Dropper cap fits for thread dia. 18 mm		0300-3018	N3
Dropper cap fits for thread dia. 25 mm		0300-3025	N3

Narrow-neck reagent bottle



Universal bottles for storing, packing, taking samples, sampling and transport of liquids.

- HDPE, transparent
- PP tamper-evident screw cap, blue

Narrow-neck reagent bottle

Capacity ml	L x W mm	Height capless mm	Thread dia. mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.
100	45 x 45	108	32	17		100		0340-0100 N3
250	60 x 60	136	32	17		50		0340-0250 N3
500	74 x 74	170	32	17		50		0340-0500 N3
1000	85 x 85	221	32	17		72		0340-1000 N3
2500	122 x 122	272	45	28		27		0340-2500 N3

Wide-neck reagent bottle



The reagent bottle is particularly suited for storing, packing, taking samples, sampling and transport of granules and paste-like media.

- HDPE, transparent
- PP tamper-evident screw cap, blue
- UN approval for bulk goods (see page 269)

Wide-neck reagent bottle

Capacity ml	L x W mm	Height capless mm	Thread dia. mm	Neck dia. inside mm	UN approval	Price €	Quantity scale from	Scale price €	Item no.
50	36 x 36	89	28	21	no		700		0342-0050 N3
100	45 x 45	101	32	24	UN3H1/X0.3		582		0342-0100 N3
250	59 x 59	130	45	35	UN3H1/X0.8		273		0342-0250 N3
500	75 x 75	162	54	43	UN3H1/X1.5		124		0342-0500 N3
750	83 x 83	188	54	43	UN3H1/X2.5		90		0342-0750 N3
1000	85 x 85	214	54	43	UN3H1/X2.5		81		0342-1000 N3
1500	108 x 108	209	80	69	UN3H1/X4		42		0342-1500 N3
2500	122 x 122	260	80	69	UN3H1/X5		27		0342-2500 N3
4000	142 x 142	294	80	69	UN3H1/X5.5		15		0342-4000 N3

Wide-necked bottle | easy to fill and clean, compact form, unbreakable





Sturdy bottle with wide-necked opening and screw cap. Particularly suitable for viscous media, powder, granules and pastes. Leakproof.

- LDPE or PP
- LDPE coloured brown for light-sensitive liquids

Individual customer-specific printing possible - please ask us!

Wide-necked bottle

Material	Capacity ml	Colour	Dia. mm	Height mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.
LDPE	50	transparent	40	80	24	32		100		0318-0050 N3
LDPE	100	transparent	48	95	24	32		100		0318-0100 N3
LDPE	200	transparent	59	114	30	40		384		0318-0200 N3
LDPE	250	transparent	64	127	31	40		50		0318-0250 N3
LDPE	500	transparent	77	158	39	50		50		0318-0500 N3
LDPE	750	transparent	88	176	39	50		98		0318-0750 N3
LDPE	1000	transparent	94	207	50	65		50		0318-1000 N3
LDPE	1500	transparent	108	232	50	65		46		0318-1500 N3
LDPE	2000	transparent	120	252	50	65		50		0318-2000 N3
PP	100	transparent	48	95	24	32		100		0319-0100 N3
PP	250	transparent	64	127	31	40		50		0319-0250 N3
PP	500	transparent	77	158	39	50		100		0319-0500 N3
PP	1000	transparent	94	207	50	65		50		0319-1000 N3
LDPE	250	brown	64	127	31	40		50		0320-0250 N3
LDPE	500	brown	77	158	39	50		50		0320-0500 N3
LDPE	1000	brown	94	207	50	65		50		0320-1000 N3

Wide-necked square bottle



Cost-effective square bottle for packing and storing of liquids, powders and granules.

- HDPE, transparent, shatterproof
- Special labelling area for labels or manual lettering
- Leakproof

Wide-necked square bottle

Capacity ml	L x W mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
250	59 x 59	110.5	30	50		100		0338-0250	N3
500	68 x 68	146	30	50		100		0338-0500	N3
1000	90 x 90	170	50	65		50		0338-1000	N3

Multi-Bottles with graduation



Bottle for universal use. The medium-sized bottle opening ensures good filling and safe pouring.

- Bottle HDPE, lid PP
- Embossed contents scale
- Absolutely leakproof
- With screw cap

Bulk prices

Multi-Bottles with graduation

Capacity ml	Container dia. mm	Height capless mm	Neck dia. inside mm	Division ml	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
100	46	99.5	19	25	25		20		0301-0100	N3
250	61	129.5	24	50	28		10		0301-0250	N3
500	76	161	30.8	50	40		10		0301-0500	N3
1000	97	203	39.2	100	50		4		0301-1000	N3

Sealable wide-necked bottles







Sturdy bottle with sealing eyelets on the screw cap and the bottle. Particularly suited for sending samples. Seals provide increased sample authenticity. Easy to fill due to wide opening. Particularly suitable for viscous media, powder, granules and pastes.

- HDPE
- Leakproof
- Easy to fill due to wide-necked opening


Sealable wide-necked bottles

Capacity ml	Container dia. mm	Neck dia. inside mm	Height with cap mm	Price €	Quantity scale from	Scale price €	Item no.	
250	64	26.5	118		50		0322-0250	N3
500	78	40	147		50		0322-0500	N3
1000	100	50	186		50		0322-1000	N3
2000	120	50	230		50		0322-2000	N3

Accessories

Description	Price €	Item no.	
SampleSafe seal, red, 100 pieces		5399-2902	N4

Spray bottle with pump vaporizer





LaboPlast

Spray bottles with pump vapouriser and safety cap for the pump. Liquids can be finely vapourised at the push of a button in an environmentally compatible way. The transparent bottles make identification of contents easy. The protective cover prevents unintentional pressing of the pump and acts as a dust cover.

- Bottle made of HDPE, spray pump made of PP and stainless steel

Spray bottle with pump vaporizer

Capacity ml	Dia. mm	Height with cap mm	Price €	Quantity scale from	Scale price €	Item no.
20	28	85		100		0309-3002 N3
50	35	110		100		0309-3005 N3
100	37	161		100		0309-3100 N3
250	50	222		100		0309-3025 N3

Spray bottle



Pressureless spray bottle made of LDPE with long-life hand pump. The functioning of each individual sprayer is tested during manufacture.
Areas of application: Industry, laboratories, production, leisure, household. Suitable for agricultural pesticides and detergents.

- Media-contacting parts PP, PE and stainless steel
- Nozzle 0.6 mm diameter, nozzle steplessly adjustable
- Spray volume per stroke: 1.2 ml ± 0.1

LaboPlast

Spray bottle

Capacity ml	Height with cap mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.
250	220	18		180		0309-1025 N4
500	240	18		80		0309-1050 N4
1000	295	17.5		50		0309-1100 N4

Spray bottle Turn ‘n’ Spray | Spray bottle with overhead valve





- Uninterrupted spraying, even over your head
- Parts in contact with media made of PP, PE, silicon and stainless steel
- Nozzle diameter 0.6 mm, nozzle is continuously adjustable
- Amount sprayed per stroke: 1.2 ml ± 0.1

Spraying over your head is not possible with standard spray bottles. As soon as the suction tube is no longer immersed in liquid, you can’t spray any more.
No so with our overhead spray bottle Turn’n’Spray. It features a special technology that allows you to spray from any position. That means you can spray practically and conveniently from any direction and reach hard-to-access places easily and without straining.
The high-quality, refillable quality spray bottle can be used again and again.

Spray bottle Turn ‘n’ Spray

Capacity ml	Height with cap mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.
250	220	18		180		0309-0002 N4
500	240	18		80		0309-0005 N4
1000	295	17.5		50		0309-0010 N4

Pressure sprayer



Pressure sprayer, robust industrial quality for disinfection, cleaning, plant protection, thin layer chromatography, technical and sanitary use. The hand pump is used to produce a pressure of up to 4 bar in the container. Spraying can then take place by pressing a button.

- Media-contacting parts PP, PE, glass-fibre-reinforced PA, V2A, FKM
- Bottle made of HDPE
- Nozzle diameter 0.6 mm
- Adjustable spray jet
- Graduation up to 1500 ml
- Temperature max. 30 °C
- Spraying lance extension 30 cm available as accessory

Pressure sprayer

Capacity ml	Container dia. mm	Height with cap mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.	
1500	142	315	45		10		0309-0100	N4

Accessories

Description	Price €	Item no.	
Spraying lance extension for pressure sprayer, 30 cm		0309-0101	N4

Wash bottle



The swan-neck jet produces a fine, straight spray with no lateral spatter. Dangerous lateral spatter does not occur.

- LDPE
- Transparent
- Immersion tube can be pulled up

LaboPlast

Wash bottle

Capacity ml	Dia. mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
100	45	200		100		0310-0010	N3
250	61	225		100		0310-0025	N3
500	75	270		100		0310-0050	N3
1000	98	315		50		0310-0100	N3

Wide-necked wash bottle



Wash bottle with medium overprint as well as hazard and safety instructions (health statements/physical statements) based on the GHS system with substance name in German, English, French and Spanish, GHS pictograms and CAS number.


- Safe filling due to wide-necked opening 39 mm neck diameter
- LDPE yellow transparent with light protection
- Contents 500 ml

LaboPlast

Wide-necked wash bottle						
Print	Colour	Price €	Quantity scale from	Scale price €	Item no.	
without label - neutral	yellow		100		0310-3050	N3
Distilled water	transparent		100		0310-3051	N3
Methanol	green		100		0310-3052	N3
Ethanol	orange		100		0310-3053	N3
Isopropanol	blue		100		0310-3054	N3
Acetone	red		100		0310-3055	N3

Safety wash bottle





The vapour-permeable membrane reliably prevents leaking of the contents.

If the room temperature increases, or the sun shines into the lab, conventional wash bottles will drip because of the increased inner pressure. The venting valve ensures that the excess pressure can escape and dripping is prevented. With leakproof venting valve, contents cannot leak when bottle has fallen over as with wash bottles with ball valve or venting hole.

- Safe filling due to wide-necked opening 39 mm neck diameter
- Permanent print according to GHS with substance name in German, English, French and Spanish, GHS pictograms and CAS number
- LDPE yellow transparent with light protection
- Contents 500 ml

LaboPlast

Safety wash bottle						
Print	Colour	Price €	Quantity scale from	Scale price €	Item no.	
without label - neutral	yellow		12		0310-2056	N3
Distilled water	transparent		12		0310-2050	N3
Methanol	green		12		0310-2052	N3
Ethanol	orange		12		0310-2054	N3
Isopropanol	blue		12		0310-2055	N3
Acetone	red		12		0310-2051	N3

Dropping bottles | for controlled dosing drop by drop



Dropping bottles are particularly suitable for liquids which must be dispensed in small amounts, drop by drop or as a fine jet.

- Cut off dropper tip before first use
- Can be reclosed using the captive cap provided
- LDPE
- Captive cap and dropper cap PE

LaboPlast

Dropping bottles

Capacity ml	Container dia. mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.
10	26	76	9	14		100		0306-0010 N3
20	32	88	9	14		100		0306-0020 N3
30	35	96	9	14		50		0306-0030 N3
50	40	116	13	18		50		0306-0050 N3
100	45	136	13	18		25		0306-0100 N3
250	61	170	18	25		25		0306-0250 N3
500	75	204	18	25		20		0306-0500 N3
1000	95	248	21	28		10		0306-1000 N3

Drop-Boy



Bottle with drop dispenser and separate screw cap.

- Cut off dropper tip before first use
- Reclosable
- Precise drop dosing

Drop-Boy

Capacity ml	Dia. mm	Height capless mm	Price €	Quantity scale from	Scale price €	Item no.
10	24	50		100		0307-0010 N3
25	31.3	64.5		100		0307-0025 N3
50	40	77		100		0307-0050 N3

UN-Bottle



This bottle reduces packaging costs! UN tested and approved. Bottle made of HDPE, leakproof, complete with PR tamper-evident closure with PTFE seal. The bottles are excellent for packaging purposes, sampling, sending samples, storage and shipping via the postal service.

- Absolutely leakproof
- With UN approval for liquids (see page 269)
- Transparent
- Bell-shaped caps with or without degassing on request

UN-Bottle

Capacity ml	Height capless mm	Neck dia. inside mm	UN approval	Price €	Quantity scale from	Scale price €	Item no.	
500	177	29	1H1/X1.9/250		120		0304-0500	N3
1000	212	29	3H1/X1.9/250		64		0304-1000	N3
2500	282	29	3H1/X1.9/250		26		0304-2500	N3

Bottle carrier



Bottle carrier
For 6 bottles up to max. 95 mm Ø.

- PP
- Suitable for food



Bottle carrier box
For 8 bottles up to max. 95 mm Ø. With closed floor which can be used as spill trough.

- HDPE
- Inside measurements: Length 390 mm, width 110 mm


Bottle carrier

No. of bottles	Container dia. mm	Outer dimension of rim L x W x H mm	Price €	Quantity scale from	Scale price €	Item no.	
6	95	310x210x305		10		0361-0006	N3

Bottle carrier box

No. of bottles	Container dia. mm	Outer dimension of rim L x W x H mm	Price €	Quantity scale from	Scale price €	Item no.	
8	95	420x300x210		10		0361-0008	N3

Aluminium bottle | light, unbreakable, seamless





The aluminium bottle is particularly suited for storing, transport or sending substances and reference samples which have to be stored completely unaltered. Other properties: can be completely emptied without residues, odour-neutral, impervious to light, rustproof.

- Pure aluminium (material AL 99.5)
- With tamper-evident screw cap made of PP with aluminium inner seal
- Suitable for food (exception: highly acidic and salt-containing foods)
- With UN approval (for 60–5600 ml), see page 269

Aluminium bottle

Capacity ml	Dia. mm	Height with cap mm	UN approval	Price €	Quantity scale from	Scale price €	Item no.
38	35	70	no		387		0327-0038 N3
60	40	88	1B1/X/250		304		0327-0060 N3
120	57	97	1B1/X/250		143		0327-0120 N3
300	57	160	1B1/X/250		149		0327-0300 N3
600	74	192	1B1/X/250		84		0327-0600 N3
1200	88	248	1B1/X/250		60		0327-1200 N3
3000	161	215	1B1/X/250		168		0327-3000 N3
5600	161	340	1B1/X/250		168		0327-5600 N3



The UN code is printed onto the bottles in black.



With tamper-evident screw cap made of PP with aluminium inner seal.



Aluminium transport bottles with screw cap close safely and tightly and are approved for packaging type I (X). The UN code is printed onto the bottles in black.

- Pure aluminium (material AL 99.5)
- UN approved (see page 269)
- With eyelets, sealable

Transport bottle aluminium UN

Capacity ml	Dia. mm	Height mm	UN approval	Price €	Quantity scale from	Scale price €	Item no.
125	50	105	1B1/X/250		10		0329-0125 N3
250	65	140	1B1/X/250		50		0329-0250 N3
500	80	145	1B1/X/250		30		0329-0500 N3
1250	100	205	1B1/X/250		45		0329-1250 N3

Tubes with screw cap



Shapely packing container with screw cap. Emptying without residues because there are no crevices, cost-effective.

- PP, white
- Sterilisable

Wide-necked boxes, round



Wide-necked containers with very good shatter resistance for convenient filling. Particularly suited for storing or shipping powders and granules.

- HDPE, transparent, shatterproof
- Screw lid made of PP with seal made of PE

Bulk prices

Tubes with screw cap

Capacity ml	Dia. mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.
12	25	31		100		6214-0012 N3
25	31	40		100		6214-0025 N3
60	41	56		100		6214-0060 N3
90	52	52		100		6214-0090 N3
125	52	67		100		6214-0125 N3
250	68	79		100		6214-0250 N3
625	90	113		70		6214-0600 N3
1000	102	129		80		6214-1000 N3

Wide-necked boxes, round

Capacity ml	Dia. mm	Neck dia. inside mm	Height capless mm	Price €	Quantity scale from	Scale price €	Item no.
100	48	31	77		100		6282-0100 N3
250	65	41	98		100		6282-0250 N3
500	80	53	122		100		6282-0500 N3
750	94	53	130		108		6282-0750 N3
1000	100	67	153		75		6282-1000 N3
1500	114	67	173		48		6282-1500 N3
2000	126	67	192		40		6282-2000 N3

Wide-necked boxes, square shaped | clear or brown material





The wide opening of the container enables easy filling with liquid, powder or granules such as tablets, capsules, dyes, shipments of samples, specimens etc. Also suitable for gift packages and hobby articles. For light-sensitive material, brown transparent wide-necked containers are particularly suited.

- PETG crystal-clear
- PVC highly transparent or brown transparent
- Screw lid PP, black, seal made of PE

Wide-necked boxes, square shaped

Material	Capacity ml	Length mm	Width mm	Height with cap mm	Neck dia. inside mm	Price €	Quantity scale from	Scale price €	Item no.
PETG crystal clear	50	38	38	63	23		162		0357-0050 N3
PETG crystal clear	100	46	46	75	31		80		0357-0100 N3
PETG crystal clear	200	59	59	91	41		100		0357-0200 N3
PETG crystal clear	500	80	80	111	54		100		0357-0500 N3
PETG crystal clear	1000	95	95	152	69		96		0357-1000 N3
PVC highly transparent	50	38	38	63	23		162		0355-0050 N3
PVC highly transparent	100	46	46	75	31		80		0355-0100 N3
PVC highly transparent	200	59	59	91	41		100		0355-0200 N3
PVC highly transparent	300	66	66	98	41		100		0355-0300 N3
PVC highly transparent	500	80	80	111	54		100		0355-0500 N3
PVC highly transparent	750	92	92	120	69		114		0355-0750 N3
PVC highly transparent	1000	95	95	152	69		96		0355-1000 N3
PVC highly transparent	2000	120	120	190	69		44		0355-2000 N3
PVC amber transparent	100	46	46	75	31		80		0356-0100 N3
PVC amber transparent	200	59	59	91	41		100		0356-0200 N3

All-purpose tubes with screw cap



Contents can be recognised easily due to the transparent sides. Emptying without residues.

- PP, transparent, screw lid PE
- No crevices

All-purpose tubes with screw cap; transparent

Capacity ml	Container dia. mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
250	104	50		25		6202-0250	N3
500	104	100		14		6202-0500	N3
1000	126	126		12		6202-1000	N3

All-purpose boxes, square



Box for storing, packing, taking samples, sampling and transport of pastes, granules, powders, soil samples and frozen goods (-40 °C up to briefly +95 °C). Ideal use of space, can be stacked. Complete with press-in lid.

- PE, transparent
- Can be deep frozen
- No crevices

All-purpose boxes, square shaped

Capacity ml	L x W x H mm	Price €	Quantity scale from	Scale price €	Item no.	
500	103 x 103 x 64		20		6250-0500	N3
1000	103 x 103 x 127		16		6250-1000	N3
1000	208 x 103 x 64		18		6250-1001	N3
1500	208 x 103 x 94		20		6250-1500	N3
2000	208 x 208 x 64		20		6250-2000	N3
2000	255 x 103 x 94		20		6250-2300	N3
3200	208 x 208 x 94		14		6250-3300	N3

All-purpose boxes



All-purpose boxes, round

Capacity ml	Outer dia. bottom mm	Outer dia. top mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
250	60	78	90		60		6210-0250	N3
500	75	95	127		24		6210-0500	N3
1000	118	136	110		20		6210-1001	N3
2000	148	169	143		25		6210-2000	N3

Sample boxes, aseptic







Sample box with tightly closing sealing cover. The lid can be opened with one hand and its hinge angle is 90°. Liquid/condensate in the lid runs into an absorbing groove integrated in the lid, nothing is lost, nothing can drip. The box is completely tight, has no crevices, is therefore easy to fill and can be emptied without residues. The sample box is produced under aseptic conditions. It is autoclavable and suitable for deep freezing.

- PP, transparent
- Aseptic
- With area for labelling and graduation

Bulk prices

Sample boxes, aseptic

Capacity ml	Height mm	Dia. mm	Model	Price €	Quantity scale from	Scale price €	Item no.	
50	78	31	with labelling area and graduation		650		6206-0060	N3
90	80	43	with labelling area and graduation at 80 ml		350		6206-0100	N3
300	88	75	graduated		240		6206-0300	N3

Sample tubes aseptic, with seal



The sample boxes are manufactured in the cleanroom and sealed in the injection moulding machine at over 100 °C, so that the insides of the boxes are virtually sterile. This method makes the boxes absolutely leakproof, as the cap and box adapt ideally to each other during the cooling process. No crevices, therefore easy to fill and can be emptied without residues. With tamper-evident seal which guarantees virtual sterility until the first opening. After opening, it can be sealed again with a practical seal.

- PP, transparent
- Aseptic
- With graduation
- Sealable

Sample tubes aseptic, with seal

Capacity ml	Height mm	Dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
55	84	29		90		6207-0055	N3
340	119	63		90		6207-0340	N3

Sample container, round



Sample box with screw lid. Ideal for taking samples and transport. Available in two material qualities: highly transparent PP for powder, pastes, liquids, or PFA for super pure collection of samples. Emptying without residues because there are no crevices.

- PP or PFA autoclavable
- Thick sides
- Sample container PFA not suitable for food

Sample container, round

Material	Capacity ml	Inner dia. mm	Outer dia. mm	Height with cap mm	Price €	Quantity scale from	Scale price €	Item no.	
PP	60	35	50	93		10		6203-0060	N3
PP	180	50	65	119		10		6203-0180	N3
PFA	90	50	65	65		10		6205-0090	N3
PFA	180	50	65	119		10		6205-0180	N3

Sample beaker with cover



Sample beaker with liquid-tight red screw-on cover. The sample beakers are aseptic and suitable for use with food. The beakers can be autoclaved without the cover.

- PP beaker, LDPE cover
- Temperature resistance of beaker -20 °C/+120 °C, cover -80 °C/+75 °C
- Aseptic

DispoDipper





Sampling tube manufactured in a clean room specially designed for disposable use. Due to the blue colour, it is also perfect for use in the food, animal feed and pharmaceutical industries. DispoDipper is completely manufactured in one piece and can be used for sampling and subsequent transport of the sample in one device.

- PP, blue transparent
- Individually packaged and optionally sterilised by gamma rays
- Production and packaging according to clean room class 7
- Complies with EU food requirements and FDA regulations
- Scaled sampling tube up to 50 ml / 1.5 oz.
- Immersion depth 270 mm

LaboPlast SteriPlast

Sample beaker with cover

Capacity ml	Dia. mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
125	52	74		380		7103-0125	N3

DispoDipper

Capacity ml	Dia. mm	Height mm	Total length mm	Model	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
50	40	98	417	standard	20		7		5393-7000	N4
50	40	98	417	sterilized	20		7		5393-7001	N4

Compact jerrycan | for more content on less shelf width



The compact jerrycan is ideal as a space-saving canister with its width of 6.5 cm or 12.5 cm.
Different liquids can be stored and decanted in a small space. The compact jerrycan can also be stored lying on its side e. g. in a refrigerator. It is very stable due to ribbed structure and thick sides.
Using the support stand (only for 5 l canister) for two or three canisters, any number of canisters can safely be placed in the open. The compact jerrycan can also be used as a settling canister.

Compact jerrycan | for more content on less shelf width









- With or without threaded connector 3/4"
- PP, transparent
- With tamper-proof closure
- Canister without stopcock autoclavable
- With moulded, indestructible litre graduation
- Canister support made of aluminium for secure stand as accessory (only for 5 l canister)
- Spigot StopCock with rotatable discharge angle and Compact spigot, which fits canisters equipped with glands with threads 3/4", available as accessory

Compact jerrycan

Model	Capacity l	Width mm	Depth mm	Height mm	Opening dia. mm	Price €	Quantity scale from	Scale price €	Item no.
without threaded connector	5	65	335	335	41		10		0435-1050 N4
with threaded connector 3/4"	5	65	335	335	41		10		0435-2050 N4
without threaded connector	10	125	335	335	41		10		0435-1100 N4
with threaded connector 3/4"	10	125	335	335	41		10		0435-2100 N4

Accessories

Description	Price €	Item no.
① Support stand for two compact jerrycans (only for 5 l canister)		0435-9102 N4
① Support stand for three compact jerrycans (only for 5 l canister)		0435-9103 N4
② Venting valve in screw cap		0435-9200 N4
③ StopCock angled (see page 44)		0525-1000 N4
④ Compact spigot (see page 45)		0530-1000 N4
Closure cap 3/4"		0600-0090 N4
⑤ Set consisting of one compact jerrycan 5 l and one spigot Compact		0435-3000 N4
⑥ Set consisting of one compact jerrycan 5 l and one spigot StopCock		0435-0001 N4
⑦ Set consisting of one compact jerrycan 10 l and one spigot Compact		0435-3010 N4
⑧ Set consisting of one compact jerrycan 10 l and one spigot StopCock		0435-0010 N4

Jerrycans





Screw lid with drain cock on jerrycan, accessory

Canister for storage or transport. Canister with UN transport approval and tamper-evident seal.

- HDPE
- UN approval 3H1/X1.9/250 (see page 269)
- Screw lid with drain cock available as accessory

Bulk prices

Jerrycans

Capacity l	L x W x H mm	Weight g	Neck dia. inside mm	Thread	Colour	Price €	Quantity scale from	Scale price €	Item no.
5	190 x 145 x 251	250	40.5	DIN50	transparent		16		1427-0005 N3
10	230 x 190 x 310	450	40.5	DIN50	transparent		10		1427-0010 N3
20	290 x 246 x 385	950	48.5	DIN60	blue		6		1427-0020 N3
30	320 x 284 x 435	1250	48.5	DIN60	blue		6		1427-0030 N3

Accessories

Description	Price €	Item no.
Screw closure DIN50 with spigot for jerrycans		1427-9102 N3
Screw closure DIN60 with spigot for jerrycans		1427-9100 N3
Original screw cap DIN50		1427-9105 N3
Original screw cap DIN60		1427-9106 N3

Electrically conductive canister with UN approval

Canister made of electrically conductive HDPE for preventing ignition sources. The canisters are capable of diverting a charge (surface resistance $<10^9$ ohms) and can be earthed, which makes them suitable for use with flammable liquids.

- HDPE, capable of diverting a charge
- With carrying handle and screw cap
- With UN approval for transport (see page 269)
- Stackable

⚠ Note the "List of Chemical Resistance"!







Electrically conductive canister with UN approval

Capacity l	L x W x H mm	UN approval	Thread	Weight g	Price €	Quantity scale from	Scale price €	Item no.
10	300 x 200 x 258	3H1/Y1.5/200	DIN 61	770		16		1428-0010 N4
20	300 x 200 x 495	3H1/Y1.5/150	DIN 61	1320		10		1428-0020 N4
30	364 x 240 x 455	3H1/Y1.5/150	DIN 61	1350		5		1428-0030 N4

Wide-necked jerrycans with/without threaded connector

Wide-necked canister with extremely wide filling opening (88 mm), easy to clean on the inside, which is particularly important when storing sensitive substances. Carrying handles and screw lid guarantee safe and convenient transport. The wide-necked canister can also be used as a settling canister. Via a threaded connector just above the base, a spigot or other accessories with an inner thread of 3/4" can be attached.

- HDPE, transparent
- UV stabilised
- Self-adhesive blue litre graduation and spigots available as accessory







Wide-necked jerrycans with/without threaded connector

Model	Capacity l	Width mm	Depth mm	Height mm	Ventilation dia. mm	Price €	Quantity scale from	Scale price €	Item no.
without threaded connector	5	130	290	280	18		6		0431-0005 N3
without threaded connector	10	137	310	400	18		5		0431-0010 N3
without threaded connector	20	165	365	485	-		3		0431-0020 N3
without threaded connector	30	240	365	485	-		2		0431-0030 N3
with threaded connector 3/4"	5	130	290	280	18		6		0431-1005 N3
with threaded connector 3/4"	10	137	310	400	18		5		0431-1010 N3
with threaded connector 3/4"	20	165	365	485	-		3		0431-1020 N3
with threaded connector 3/4"	30	240	365	485	-		2		0431-1030 N3

Accessories

Description	Price €	Item no.
Graduation for wide-necked jerrycans 5 l		0431-9051 N3
Graduation for wide-necked jerrycans 10 l		0431-9011 N3
Graduation for wide-necked jerrycans 20 l		0431-9021 N3
Graduation for wide-necked jerrycans 30 l		0431-9031 N3
StopCock angled (see page 44)		0525-1000 N4
Compact spigot (see page 45)		0530-1000 N4
Closure cap 3/4"		0600-0090 N4

Storage bottles, with/without threaded connector

Robust storage bottle with particularly thick sides with sturdy carrying handle for safe transport and storage. Blue litre graduation printed on. The venting valve (accessory) lets air into the container in order to prevent a vacuum in the bottle. The canister can also be used as a settling canister. Via a threaded connector, a spigot or other accessories with an inner thread/cap nut of 3/4" can be attached.

- HDPE
- Bottle opening diameter 50 mm







Storage bottles, with/without threaded connector

Model	Capacity l	Container dia. mm	Height capless mm	Handles	Price €	Quantity scale from	Scale price €	Item no.
without threaded connector	5	167	330	1		11		0401-0005 N3
without threaded connector	10	206	427	1		6		0401-0010 N3
without threaded connector	25	278	580	2		5		0401-0025 N3
without threaded connector	60	346	830	2		4		0401-0060 N3
with threaded connector	5	167	330	1		11		0402-0005 N3
with threaded connector	10	206	427	1		6		0402-0010 N3
with threaded connector	25	278	580	2		5		0402-0025 N3
with threaded connector	60	346	830	2		4		0402-0060 N3

Accessories

Description	Price €	Item no.
Ventilation seal for storage bottles		0401-9100 N3
① StopCock angled (see page 44)		0525-1000 N4
② Compact spigot (see page 45)		0530-1000 N4

Wide-mouth containers with handle



Container with 2.3 l and 4.4 l contents and particularly wide opening. Therefore ideal for filling with bulk goods of all kinds. Products such as spices, tea, granulates can be packaged, stored or displayed as samples in these containers. Space-saving square shape.

- Safe handling with integral handle
- Transparent polyethylene provides good view of contents.
- Leakproof wide-necked lid
- 0317-2300 with tamper-evident screw cap made of blue PP, 0317-4400 with screw cap made of LDPE

Wide-mouth containers with handle

Capacity ml	L x W mm	Height with cap mm	Neck dia. inside mm	Thread dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
2300	135 x 135	192	68.5	80		36		0317-2300	N3
4400	142 x 142	263	109	120		52		0317-4400	N3

Universal storage containers



Practical storage containers with closed gripping edges for safe transport of pourable or liquid contents. Rounded inner edges and smooth inner surfaces for easy cleaning enable use in the foodstuffs/hygiene sector.

- HDPE, white
- Can be stacked

⚠ Attention! Make sure to pay attention to the "List of Chemical Resistance"!

Universal storage containers

Capacity l	Outer dimension of rim L x W x H mm	Inner dimension L x W x H mm	Price €	Quantity scale from	Scale price €	Item no.	
29	660 x 450 x 130	600 x 400 x 120		10		3418-0129	N2
48	660 x 450 x 220	600 x 400 x 210		10		3418-0148	N2
94	660 x 450 x 410	600 x 400 x 400		10		3418-0194	N2

Accessories

Description	Price €	Item no.	
Lid for universal storage containers, PS white		3418-0001	N2

Storage and stacking containers



Stable storage containers made of durable PP with external reinforcement ribs.

- Temperature resistant from -40 °C to +90 °C
- According to Euro standards
- Lid and transport rollers available as accessory



Transport roller



Lids for storage and stacking containers

Storage and stacking containers

Capacity l	Outer dimension of rim L x W x H mm	Inner dimension L x W x H mm	Price €	Quantity scale from	Scale price €	Item no.
11	400 x 300 x 120	367 x 267 x 116		10		3414-0111 N2
17	400 x 300 x 180	367 x 267 x 176		10		3414-0117 N2
23	400 x 300 x 220	367 x 267 x 216		10		3414-0120 N2
23	600 x 400 x 120	567 x 367 x 114		10		3414-0123 N2
38	600 x 400 x 180	567 x 367 x 176		10		3414-0138 N2
45	600 x 400 x 220	567 x 367 x 216		10		3414-0148 N2
64	600 x 400 x 320	567 x 367 x 316		10		3414-0164 N2
80	600 x 400 x 420	567 x 367 x 416		10		3414-0181 N2

Accessories

Description	Price €	Item no.
Lid for storage/stacking containers 400 x 300 mm		3414-0001 N2
Lid for storage/stacking containers 600 x 400 mm		3414-0002 N2

Transport trolley for containers

Model	Price €	Quantity scale from	Scale price €	Item no.
for 600 x 400 mm		10		3414-9900 N2

Tubs, round



Tub for every purpose, with strong edge and practical carrying handles for safe transport. A lid can be delivered for every tub.

- Tub made of HDPE
- Temperature resistant from -30 °C to +70 °C
- Suitable for food
- High degree of chemical resistance

⚠ Attention! Make sure to pay attention to the "List of Chemical Resistance"!

Mixing/dosing container



The containers close completely tight. They have thick sides, they are seamless, impact resistant and stressfree sintered. The wide-necked opening enables safe filling.

- LDPE
- With 3/4" inner thread at the container base
- Opening Ø 150 mm
- Suitable for food

Tubs, round

Capacity l	Dia. top/below x H mm	Price €	Item no.	
40	370/330 x 440		3482-0040	N3
50	370/330 x 555		3482-0050	N3
75	460/400 x 535		3482-0075	N3
110	460/420 x 710		3482-0110	N3

Accessories

Description	Price €	Item no.	
Lids for tubs, round, for 370 mm dia.		3482-0401	N3
Lids for tubs, round, for 460 mm dia.		3482-0751	N3

Mixing/dosing container

Capacity l	Dia. mm	Height mm	Division l	Price €	Item no.	
60	450	600	10		3460-0060	N2
100	465	850	10		3460-0100	N2
200	560	1020	20		3460-0200	N2
500	815	1140	50		3460-0500	N2

Disposal kegs, wide-mouth containers



Wide-necked container for disposal or storage of powders, granulates and liquids. Can be closed tightly and sealed. The containers have an UN transport approval for solids.

- HDPE
- Temperature resistant from -20 °C to +80 °C
- UN transport approval for solids 1H2/X (see page 269)

Disposal kegs, wide-mouth containers

Capacity l	Container dia. mm	Height with cap mm	Opening dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
3,6	198	171	136		10		1105-0003	N3
6	198	263	136		10		1105-0006	N3
10	274	238	203		10		1105-0010	N3
26	316	425	203		10		1105-0026	N3
34	338	449	295		10		1105-0033	N3

Accessories

Description	Price €	Item no.	
SampleSafe seal, red, 100 pieces		5399-2902	N4

Wide-mouth containers, HDPE



Storage or transport container made of HDPE. For liquids, powders or granules. Can be closed tightly.

⚠ Keep the container upright when transporting liquids!

Wide-mouth containers, HDPE

Capacity l	Container dia. mm	Height with cap mm	Opening dia. mm	Price €	Quantity scale from	Scale price €	Item no.	
5	192	250	110		14		1406-1005	N3
10	232	325	110		6		1406-1010	N3

Containers, rectangular and round shaped



⚠ Depending on temperature and the expected contents, the rectangular containers may have to be strengthened with a frame. Kindly you ask us.

LaboPlast

Container for various applications. Produced seamlessly and stress-relieved with the LaboPlast® sintering method. The wall thickness increases downwardly. This guarantees optimum ratio of wall thickness to filling load. The containers are also colored black for UV-protection and are available in other wall thicknesses.
Custom-made, exchange excluded.

- L-LDPE (linear LDPE), transparent
- Rupture-proof

We will mount the container connector directly into the container for you. The thread connectors are possible anywhere. Please state the exact position on the container, if applicable.

Container thread

Model	Thread	Price €	Item no.
Container thread	3/4"		3440-0034 N2

Accessories

Description	Price €	Item no.
StopCock angled (see page 44)		0525-1000 N4
Compact spigot (see page 45)		0530-1000 N4

Containers, rectangular

Capacity l	Inner dimension L x W x H mm	Wall thickness mm	Price €	Item no.	
13	220 x 220 x 410	5		3440-0135	N2
23	540 x 370 x 130	5		3440-0235	N2
40	390 x 190 x 620	5		3440-0405	N2
45	600 x 600 x 165	5		3440-0455	N2
48	400 x 400 x 300	5		3440-0485	N2
50	1000 x 320 x 160	5		3440-0505	N2
56	315 x 315 x 635	5		3440-0565	N2
70	460 x 310 x 570	5		3440-0705	N2
85	600 x 340 x 480	5		3440-0855	N2
90	600 x 400 x 410	6		3440-0906	N2
100	780 x 460 x 320	6		3440-1006	N2
115	500 x 400 x 600	6		3440-1156	N2
118	600 x 250 x 820	6		3440-1186	N2
130	710 x 330 x 560	6		3440-1306	N2
150	600 x 400 x 630	6		3440-1506	N2
155	600 x 500 x 560	6		3440-1556	N2
165	780 x 500 x 420	6		3440-1656	N2
190	600 x 400 x 800	6		3440-1906	N2
250	1000 x 600 x 450	6		3440-2506	N2
260	800 x 600 x 600	7		3440-2607	N2
285	600 x 600 x 820	7		3440-2857	N2
320	1000 x 400 x 800	7		3440-3207	N2
400	860 x 560 x 840	7		3440-4007	N2

Lid for containers, rectangular

Model	Price €	Item no.	
for 13 l container		3440-0131	N2
for 23 l container		3440-0231	N2
for 40 l container		3440-0401	N2
for 45 l container		3440-0451	N2
for 48 l container		3440-0481	N2
for 50 l container		3440-0501	N2
for 56 l container		3440-0561	N2
for 70 l container		3440-0701	N2
for 85 l container		3440-0851	N2
for 90 l container		3440-0901	N2
for 100 l container		3440-1001	N2
for 118 l container		3440-1181	N2
for 130 l container		3440-1301	N2
for 150 l container		3440-1501	N2
for 155 l container		3440-1551	N2
for 165 l container		3440-1651	N2
for 190 l container		3440-1901	N2
for 250 l container		3440-2501	N2
for 260 l container		3440-2601	N2
for 285 l container		3440-2851	N2
for 320 l container		3440-3201	N2
for 400 l container		3440-4001	N2

Containers, round shaped

Capacity l	Inner dia. mm	Height mm	Wall thickness mm	Price €	Item no.	
14	245	325	5		3490-0144	N2
25	310	390	5		3490-0254	N2
30	310	430	5		3490-0305	N2
35	320	530	5		3490-0355	N2
50	410	480	5		3490-0505	N2
55	310	915	5		3490-0555	N2
60	410	530	5		3490-0605	N2
80	430	590	5		3490-0805	N2
85	390	710	5		3490-0855	N2
100	470	680	5		3490-1005	N2
110	540	540	6		3490-1106	N2
120	470	720	6		3490-1206	N2
145	600	510	6		3490-1456	N2
150	550	730	6		3490-1506	N2
160	460	1020	6		3490-1606	N2
180	540	820	6		3490-1806	N2
200	560	840	6		3490-2006	N2
205	800	440	6		3490-2056	N2
210	540	950	6		3490-2106	N2
250	550	1040	6		3490-2506	N2
290	740	770	6		3490-2906	N2
300	690	830	6		3490-3006	N2
400	760	1050	6		3490-4006	N2
850	970	1170	7		3490-8507	N2

Lid for containers, round shaped

Model	Price €	Item no.	
for 14 l container		3490-0141	N2
for 25 l container		3490-0251	N2
for 30 l container		3490-0301	N2
for 35 l container		3490-0351	N2
for 50 l container		3490-0501	N2
for 55 l container		3490-0551	N2
for 60 l container		3490-0601	N2
for 80 l container		3490-0801	N2
for 85 l container		3490-0851	N2
for 100 l container		3490-1001	N2
for 110 l container		3490-1101	N2
for 120 l container		3490-1201	N2
for 145 l container		3490-1451	N2
for 150 l container		3490-1501	N2
for 160 l container		3490-1601	N2
for 180 l container		3490-1801	N2
for 200 l container		3490-2001	N2
for 205 l container		3490-2051	N2
for 210 l container		3490-2101	N2
for 250 l container		3490-2501	N2
for 290 l container		3490-2901	N2
for 300 l container		3490-3001	N2
for 400 l container		3490-4001	N2
for 850 l container		3490-8501	N2

Safety jerrycans stainless steel



- 0.8 mm sides
- Handles made of solid stainless steel
- Closure thread 1 1/2"
- Seals made of PTFE
- Accredited for UN-transports according to ADR, RID,IMDG Code, ICAO-TI, GGVSEB, GGVSee

Safety canister for transport and storage made of stainless steel AISI 316 Ti (1.4571) for 5 l and 10 l canisters and AISI 304 (1.4301) for 20 l canisters. Protection from exploding through flame protection and excess pressure valve. The containers have a high degree of chemical resistance and they save space.
In canisters with excess pressure valve, pressure balance takes place from about 0.3 bar.
Safety canister with self-closing fine dosage spout 20 l: with separate ventilation. Canisters with a venting screw enable manual air supply when used or manual pressure balance during storage.

Safety cans stainless steel



- Stainless steel AISI 316 Ti (1.4571), 0.6 mm thickness
- Seals made of PTFE
- Handle made of powder-surfaced aluminium
- Closure thread 1 1/4"

For safe metering and handling of smaller amounts of flammable liquids. Protection from exploding through flame protection and excess pressure valve, the containers have a high degree of chemical resistance.
Can with removable flame protection in the filling connector and excess pressure valve in the self-closing fine dosage spout (pressure balance from about 0.3 bar).

Safety jerrycans stainless steel

Model	Capacity l	Purpose	L x W x H mm	UN approval	Price €	Item no.
with screw cap without pressure-relief valve	10	transport	280 x 130 x 380	3A1/X1.6/300		2601-7010 N2
with screw cap and pressure-relief valve	5	storage	200 x 130 x 310	no		2601-3005 N2
with screw cap and pressure-relief valve	10	storage	280 x 130 x 380	no		2601-3010 N2
with screw cap and pressure-relief valve	20	storage	345 x 175 x 470	no		2601-3020 N2
with fine dosage spout incl. pressure-relief valve	10	storage/transfer	280 x 130 x 380	no		2601-2010 N2
with fine dosage spout, pressure-relief valve, ventilation	20	storage/transfer	345 x 175 x 470	no		2601-2020 N2

Safety cans stainless steel

Capacity l	Dia. mm	Height mm	Purpose	UN approval	Price €	Item no.
1	111	235	storage/transfer	no		2601-1001 N2
2	131	295	storage/transfer	no		2601-1002 N2
5	165	400	storage/transfer	no		2601-1005 N2

Safety storage container stainless steel

Safety containers for transport and storage made of high-quality stainless steel AISI 304 (1.4301) / AISI 316 Ti (1.4571) with thickness of sides of 0.8 mm (10 l), 1.2 mm (25 l) and 1.5 mm (50 l). The 25 l container can be stacked.
Protection from exploding through flame protection and excess pressure valve, seals made of PTFE, the containers have a high degree of chemical resistance.
In safety containers with excess pressure valve, pressure balance takes place from about 0.3 bar.
For containers with self-closing tap, constant ventilation can be set via a plastic lever on the screw cap.
The contents display made of perfluoroalkoxy copolymers (PFA) enables a visual contents check when filling and emptying.

- Accredited for UN-transports according to ADR, RID,IMDG Code, ICAO-TI, GGVSEB, GGVSee



① ② ③ ④ ⑤

Safety storage container stainless steel

Model	Capacity l	Purpose	Height mm	Thread	Dia. mm	UN approval	Thread, outflowing side	Price €	Item no.
① with screw cap, without pressure-relief valve	25	transport/storage	530	1 1/2"	300	1A1/X1.6/400	-		2601-5125 N2
② with screw cap and pressure-relief valve	50	storage	600	G 2"	365	no	-		2601-5050 N2
③ with tab and pressure-relief valve	10	storage/transfer	330	1 1/2"	260	no	3/4"		2601-4010 N2
④ with tab and pressure-relief valve	25	storage/transfer	540	1 1/2"	315	no	3/4"		2601-4025 N2
⑤ with tab, pressure-relief valve and level indicator	10	storage/transfer	330	1 1/2"	260	no	3/4"		2601-6010 N2

Packaging bucket



Packaging bucket for transport, packaging and storing of various materials. Suitable for the food or building material industries and non-aggressive chemicals in environmental protection, hospitals and the household. With strong handle. Quality closure on the bucket for product safety.

- PP, white
- Air-tight closure
- Complete with lid
- With advertising printing in large quantities on request

Bulk prices

Packaging bucket

Capacity l	Dia. top mm	Height with cap mm	Price €	Quantity scale from	Scale price €	Item no.	
1	132	131		100		2327-0001	N3
2	175	132		80		2327-0020	N3
3	225	123		45		2327-0030	N3
5	225	194		36		2327-0050	N3
10	268	263		25		2327-0010	N3

Packaging bags with closure



Bags with the "zipper" closure for storing, collecting, protecting and shipping of specimens, samples and sales products. Contents remain dust-free. The film is waterproof and has good tear resistance.

- LDPE, crystal-clear
- 0.05 mm strong film

Bulk prices

Packaging bags with closure

L x W mm	Volume ml	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
100 x 70	80	100		10		2348-0002	N3
130 x 60	100	100		10		2348-0001	N3
120 x 80	140	100		10		2348-0003	N3
150 x 100	300	100		10		2348-0004	N3
170 x 120	400	100		10		2348-0005	N3
220 x 160	1000	100		10		2348-0006	N3
250 x 180	1500	100		10		2348-0007	N3
300 x 200	2500	100		10		2348-0008	N3
350 x 250	4500	100		10		2348-0009	N3
400 x 300	6500	100		10		2348-0010	N3

SteriBag StandUp sample bag







SteriBag StandUp is a sterile sample bag with a stand-up base. The special bottom unfolds as soon as the bag is filled. This makes SteriBag StandUp stable and able to stand on its own without a separate stand or support, leaving both your hands free to take the sample and fill the bag. The wire bracket at the closure keeps the bag open, making it easy and convenient to fill.
The sample bag is the ideal space-saving, unbreakable, inexpensive and sterile replacement for bottles and beakers.
Can be used for samples available as liquids, powders, solids and pastes; suitable for foods and can be frozen.
The stable standing SteriBag StandUp sample bag is easy and safe to transport.
SteriBag StandUp, the sample bag for easy handling of samples.

Bulk prices

SteriBag StandUp sample bag

Capacity ml	L x W mm	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
400	229 x 114	500		5		5344-0400	N3
600	229 x 152	500		5		5344-0600	N3

Sampling bag SteriBag



The sterilised sampling bag SteriBag for samples in chemistry, pharmaceuticals, medicine, veterinary science and foodstuffs etc.
The bags are made of polyethylene and are therefore exceptionally tearproof and leakproof. See page 130.

Sampling bag SteriBag

Capacity ml	Model	L x W mm	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.	
60	without labelling area	127 x 76	500		6		5344-0060	N3
60	with labelling area	127 x 76	500		6		5344-0061	N3
150	without labelling area	178 x 76	500		6		5344-0180	N3
150	with labelling area	178 x 76	500		6		5344-0181	N3
450	without labelling area	229 x 114	500		6		5344-0540	N3
450	with labelling area	229 x 114	500		6		5344-0541	N3
650	without labelling area	229 x 140	500		6		5344-0720	N3
650	with labelling area	229 x 140	500		6		5344-0721	N3
900	without labelling area	305 x 114	500		6		5344-0810	N3
900	with labelling area	305 x 114	500		6		5344-0811	N3
1650	without labelling area	305 x 178	250		12		5344-1650	N3
1650	with labelling area	305 x 178	250		12		5344-1651	N3