Nuclear Fast Red Solution

Nuclear Fast Red Solution ab246831 is a stain with histological applications.

Key facts

Form

Liquid

Storage

Shipped at conditions

Ambient - Can Ship with Ice

Appropriate short term storage conditions

Ambient

Appropriate long term storage conditions

Ambient

Notes

Nuclear Fast Red Solution ab246831 is a stain with histological applications. The reagent has improved stability over current formulations. Current formulations tend to precipitate in cold temperatures such as experienced during winter shipping. In addition, most formulations develop a small amount of precipitate over extended periods of time. This advanced formulation eliminates problems associated with exposure to cold and aging.Nuclei: RedCytoplasm: Pale PinkControl Tissue: Any well fixed paraffin embedded or frozen tissue section. Cell Smear. Cytospin.Procedure:1. Deparaffinize sections if necessary and hydrate to distilled water.2. Apply Nuclear Fast Red Solution (Enhanced Stability) and incubate for 1-5 minutes.3. Rinse in 2 changes of distilled water.4. Dehydrate through graded alcohols.5. Clear, and mount in synthetic resin.

Product promise

We are dedicated to supporting your work with high quality reagents and we are here for you every step of the way should you need us.

In the unlikely event of one of our products not working as expected, you are covered by our product promise. Full details and terms and conditions can be found here:
https://www.abcam.com/en-nl/legal/terms-and-conditions