1-800-518-0364 / 1-450-973-9420
1-450-973-6373
www.ga-international.com
www.labtag.com | www.industritag.com

PRODUCT DATA SHEET

Thermal-Transfer Xylene & Solvent Resistant Resin Ribbon- XAR Class Ribbons

Outstanding mechanical durability for most demanding applications. Prints razor sharp images on many premium synthetic materials.

Ink Color:	Black
Coating:	Coated side out (CSO)
Chemical Resistance:	Strong and genuine resistance to harsh chemicals including xylene, toluene, hexane, acetone, alcohols, thinners when used with XST-, HTT-, FTT-, AAA- class labels from GA International Inc.

Intended Applications:

Hospital / Medical Laboratories – labeling of microscope slides and tissue processing cassettes in histology and cytology labs, histological staining of slides, labeling of vials and other containers. Other applications in electronic, automotive, chemical and other industries where solvents and chemicals are used.



Storage Conditions:	Room temperature below 104ºF (40ºC), away from flame and heat sources.
RoHS Compliance:	According to the primary material manufacturer, no chemical listed in EU Directive 2002 / 95 / EC, as amended (Restriction of the use of certain Hazardous Substances in Electrical and Electronic Equipment), and / or in the recast RoHS Directive 2011/ 65 / EU is knowingly added or intentionally. This includes the restriction of Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP) and Diisobutyl phthalate (DIBP).

1-800-518-0364 / 1-450-973-9420
1-450-973-6373
www.ga-international.com
www.labtag.com | www.industritag.com

REACH Compliance	Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH). Certifying that a product is compliant with the EU REACH regulation (EC) No 1907/2006.
Product Warranty:	Refer to our *Terms and Conditions* found on www.labtag.com to know about *our 12-Month Warranty Against Manufacturing Defects*

The specifications provided herein are based on the information supplied by primary material manufacturer (s) and tests conducted by GA International Inc. The end-user shall independently determine the suitability of the labels for the purpose for which they were purchased.

USA
P.O. Box 2633
Champlain, NY
12919, United States

CANADA
3208 Jacques-Bureau Ave.
Laval, QC
H7P 0A9, Canada

NETHERLANDS
Runweg 28A
5258 BN Berlicum
Netherlands

UK
10 John Street
London WC1N 2EB
United Kingdom