

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3753	PCR PLATE 96W NON SKIRTED 0.2ML	0.2ml	Standard Profile	Non-Skirted	Natural	-
732-3754	PCR PLATE 96W NON SKIRTED 0.2ML WHITE	0.2ml	Standard Profile	Non-Skirted	White	-



SECTION B-B



Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code

2 HDPE 4 LDPE 5 PP

Description: PCR Plate


Drawing: Not to Scale


avantor™

Release: 3/25/2021

Revision: 0

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3758	PCR PLATE SEMISKIRT RAISED DECK 96W 0,2	0.2ml	Standard Profile	Semi Skirted Raised Deck	Natural	-
732-3759	PCR PLATE SEMISKIRT RAISDECK 96W 0,2 BC	0.2ml	Standard Profile	Semi Skirted Raised Deck	Natural	Yes
732-3760	PCR PLATE SEMISKIRT RAISDECK 96W 0,2 WH	0.2ml	Standard Profile	Semi Skirted Raised Deck	White	-





Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code

2 HDPE 4 LDPE 5 PP

Description: PCR Plate


Drawing: Not to Scale

avantor™


Release: 3/25/2021

Revision: 0

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3755	PCR PLATE 96W SEMI SKIRTED 0.2ML	0.2ml	Standard Profile	Semi Skirted	Natural	-
732-3756	PCR PLATE 96W SEMI SKIRTED 0.2ML BCODE	0.2ml	Standard Profile	Semi Skirted	Natural	Yes
732-3757	PCR PLATE 96W SEMI SKIRTED 0.2ML WHITE	0.2ml	Standard Profile	Semi Skirted	White	-



SECTION B-B



Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code

2 HDPE 4 LDPE 5 PP

Description: PCR Plate


Drawing: Not to Scale


avantor™

Release: 3/25/2021

Revision: 0

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3761	PCR PLATE 96W NON SKIRTED LOW P 0.1ML	0.1ml	Low profile	Non-Skirted	Natural	-
732-3762	PCR PLATE 96W NON SKIRTED LOWP 0.1ML WH	0.1ml	Low profile	Non-Skirted	White	-





Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code

2 HDPE 4 LDPE 5 PP

Description: PCR Plate


Drawing: Not to Scale

avantor™


Release: 3/25/2021

Revision: 0

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3768	PCR PLATE 96W SKIRTED 0.1ML	0.1ml	Low profile	Full Skirted	Natural	-
732-3769	PCR PLATE 96W SKIRTED 0.1ML BARCODED	0.1ml	Low profile	Full Skirted	Natural	Yes
732-3770	PCR PLATE 96W SKIRTED 0.1ML WHITE	0.1ml	Low profile	Full Skirted	White	-



SECTION A-A



SECTION B-B

Material Recycle Code

2 HDPE 4 LDPE 5 PP

Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Description: PCR Plate


Drawing: Not to Scale

avantor™


Release: 3/25/2021


Revision: 0

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3763	PCR PLATE 96W SEMISKIRT ABI FAST LOW P	0.1ml	Low profile	Semi Skirted ABI FAST	Natural	-
732-3764	PCR PLATE 96W SEMISKIRT ABI FAST LOW BC	0.1ml	Low profile	Semi Skirted ABI FAST	Natural	Yes
732-3765	PCR PLATE 96W SEMISKIRT ABI FAST LOW WH	0.1ml	Low profile	Semi Skirted ABI FAST	White	-



A-A





SECTION B-B

Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code



Description: PCR Plate

Drawing: Not to Scale


avantor™


Release: 3/25/2021

Revision: 0

Product Cat#	Description	Capacity	Profile Type	PCR Plate Type	Colour	Barcoded
732-3766	PCR PLATE 96W SEMISKIRT ROCHE 0.1ML WH	0.1ml	Low profile	Semi Skirted Roche type	White	-
732-3767	PCR PLATE SEMISKIRT ROCHE 96W 0.1 WH BC	0.1ml	Low profile	Semi Skirted Roche type	White	Yes







SECTION A-A

Material : 100 % Virgin polypropylene, US FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code



Description: PCR Plate

Drawing: Not to Scale

avantor™

Release: 3/25/2021

Revision: 0