Datenblatt

Technical data sheet

Sorte | grade: **LabSorb**

Art. Nr. | item No.: **601**

Eigenschaften | characteristics: Absorbent paper coated with polyethylene.

Anwendung | applications: Protection of work surfaces from spilled chemicals.

Prüfmerkmal property	Prüfmethode test method	Einheit unit	charakteristische Werte typical values
Flächenbezogene Masse Basis weight	DIN EN ISO 536	g/m²	140
Wasseraufnahme Water absorption	intern internal	%	150

23.11.2023 Rev.03

(*)Approximate value, not part of the product specification.

All indicated data to be understood as typical average values.
This product is manufactured in conformance with the Quality Management standard DIN EN ISO 9001.

Contact office:

Ahlstrom Germany GmbH
Niederschlag 1
09471 Bärenstein, Germany
T: +49 (0) 37347 / 83-0
F: +49 (0) 37347 / 83-64
filtrak@ahlstrom.com

www.ahlstrom.com

DISCLAIMER: The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as Ahlstrom deems appropriate. Refer to www.ahlstrom.com for contact information.
© Ahlstrom 2023