DATA SHEET Gr. 0604

●DESIGNATION

-Bench protective paper
-Material: 80% cellulose 20% PE film
-Binder: no

●APPLICATIONS

-Protection of laboratory benches, working areas cupboard or trays
-Protection against contamination, chemical pills, toxic infectious and radioactive materials
-Suitable for laying out of animal cages

●TECHNICAL SPECIFICATIONS

Characteristics	Values	Unities	Standards
Weight	125 (+/- 5%)	g/sqm	150536
Thickness	0,24	mm	
Capillary rise	50-60	mm	Klemm
Capacity of absorption	170	g/m²	