TECHNICAL DATA SHEET

General Catalogue 2021 page 191


2 – 10 ml macrotip

codes 200015, 200015R, 200015RF

code	description	characteristics	sterile	case quantity	case weight (kg)	case volume (m^3)
200015	tips in bulk	●	no	10 x 100	7.24	0.087
200015R	25 tips rack	●	no	4 x 25	2.25	0.015
200015RF	25 tips rack filter	● ● ●	STERILE R	4 x 25	2.20	0.015

● RNAsa and DNAsa free. ● Metal free. ● Pyrogen free.





USE AND DESCRIPTION

Natural colour universal macro tip with collar. Gilson® type.

For use (but not limited to):
Labsystem®, Socorex®, Finpipette®.

Rack made of polypropylene autoclavable up to 121 ºC.

Transparent lid and marble base.

Sliding lid.

DELTALAB S.L confirms that it is neither the owner of the brand Gilson® nor has any partnership with the company of this proprietary brand.

WRITTEN	VERIFIED	AUTHORISED
Joan Ibáñez Marketing Department	Anna Mir Technical, Development and Quality & Environment Director	Mónica Torras Sales & Marketing Director