HIRSCHMANN®

Data sheet

Article: Measuring cylinder 25 ml Class A
Code no.: 224 01 70

Material:
DURAN (Borosilicat glass 3.3)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 1 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 1 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Nominal volume:
25 ml

Subdivision of the scale:
0,5 ml

Tolerance:
± 0,25 ml

Hirschmann Laborgeräte
GmbH & Co.KG

Roland Benz
Quality Manager

Eberstadt, 25.03.2024