

Durchgang mit Kappe überprüfen!!
Check pass with cap!!


GL45 Edelstahl 316L/1.4404 Einsatz 2-Port


Silikon Ø40,5x6,5x3
60°± 5° Shore A
D294901-GZA-0GL45-..

Schraubverbindungskappe GL45
mit Loch
Art.-Nr. 29 227 10 15









poliert
(Ra 0,8)

B	Länge Anschluß war 70 (117)	16-MAY-2018	Weindorf
A	Rohr Rauigkeit Ra innen 1,6µm u.Ra außen 0,4µm hinzu	14-JUN-2017	Weindorf
Index	Aenderung / changes	Datum / date	Name

Weitergabe sowie Vervielfältigung dieses Dokuments, Verwertung und Mitteilung seines Inhalts sind verboten, soweit nicht ausdrücklich gestattet. Zuwiderhandlungen verpflichten zu Schadenersatz.
Alle Rechte für den Fall der Patent-, Gebrauchsmuster- oder Geschmacksmustereintragung vorbehalten.
The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited.
Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

zul. Abweichung allowed dev.	ISO 2768m-K	Werkstoff / material	Stück quantity / Gewicht weight / Pos. position
Oberfläche surface	ISO1302		
	Datum / date: Name	Benennung / title	Maßstab scale
bearb. drawn	26.09.2016 Weindorf	Stainless steel hose Connection 2 port Art.-Nr. 29 261 27 18; 27 999 00 18 (Laser) Art.-Nr. 29 261 27 01 mit Kappe	1 : 1
geprüft checked		Form-Nr. mould number	A3

DURAN WHEATON KIMBLE
DURAN Produktions GmbH & Co. KG
www.DWK-LifeSciences.com

Zeichnungsnummer / drawing number	Revision	Blatt/sheet
D362761-GZA-0GL45	B	1 /1
Ers. / repl.f:	Ers. / repl.d:	
Schutzvermerk ISO 16016 / Refer to protection notice ISO 16016		INV